



## 7" Digital Lcd Microscope Biological USB Portable Dual Lens 2.0M A33.5102

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)


### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A33.5102
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




### Product Specification

- LCD Screen: 7" IPS Screen, With Touch Button To Take Photo & Video, 2.0M, 30fps
- Head: Mono Straight Head
- Eyepiece: WF10x, WF16x
- Nosepiece: Triple
- Objective: Achromatic 4x10x40x
- Working Stage: Mechanical Stage
- Focusing: Coarse & Fine Focusing
- Highlight: **7" digital lcd microscope,**  
**usb digital lcd microscope,**  
**2.0m lcd microscope**



### More Images



## Product Description

7" LCD Biological + USB Portable Dual Lens Digital Microscope, 2.0MA33.5102

Dual Optical System: 7" LCD Digital Camera + USB Portable Digital Microscope  
Total Magnification 40x-640x Achromatic Optical System  
Double Layer Mechanical Stage 115x125mm With Disc Condenser  
Coarse & Fine Focusing Knob With Insert LED Light Source Brightness Adjustable  
With 10 Prepared Slides, With Experiment Kit, Pack In Hand Carry Case+Carton



Our Product Introduction


for more products please visit us on [opto-edu.com](http://opto-edu.com)

## Ergonomics

The best viewing angle range, sitting for a long time will not be tired and will not hurt your eyes

## 7" LCD Measurable

The screen can adjust the visual angle of 180°, real-time real-time measurement of cell size





### Objective & Stage

Triple, Achromatic 4x10x40x, Mechanical Stage



### Focusing

Coarse & Fine Focusing, Brass dovetail guide rail, high precision and durable



A33.5102 7" LCD Biological Digital Microscope, 2.0M	
LCD Screen	7" IPS Screen, With Touch Button To Take Photo & Video, 2.0M, 30fps
Head	Mono Straight Head
Eyepiece	WF10x, WF16x
Nosepiece	Triple
Objective	Achromatic 4x10x40x
Working Stage	Mechanical Stage
Focusing	Coarse & Fine Focusing
Illumination	Insert LED
TF Card Output	Support Max 32G, To Take Photo & Video
USB Output	To Win/Mac, With Professional StrangeView Measuring Software
Power Supply	DC Adapter, Can Be Powered By USB Power Bank For Outdoor Use
Package	Color Box + Foam, 1 pc/1 Box, With Software Disc, DC Adapter, Instruction Manual



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China