



Binocular Achromatic Metallurgical Optical Microscope A13.2704

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A13.2704
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Learning, Metallurgical
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Magnification: 100x - 800x
- Theory: Metallurgical Optical Microscope
- Highlight: **reflected light microscope,**
halogen lamp microscope



More Images



for more products please visit us on cnoec.com

Product Description

Binocular Achromatic Metallurgical Optical Microscope A13.2704

Specifications:

Specification	A13.2704 Inverted Metallurgical Microscope	
Optical System	UCIS infinity independent achromatic optical system	
Frame	One-piece construction from arm to base	
Head	Seidentopf trinocular head, 30° inclined, interpupillary distance 52-75mm, diopter adjustable	
Eyepiece	WF10x/20mm, plan, high eyepoint up to 21mm,	
	WF10x/18mm, high eyepoint up to 21mm,, reticule eyepiece, diopter adjustable	
Objective	Infinity LWD plan achromatic objective 10x	
	Infinity LWD plan achromatic objective 20x(s)	
	Infinity LWD plan achromatic objective 50x(s)	
	Infinity LWD plan achromatic objective 80x(s)	
Focus System	Coaxial coarse and fine focus mechanism with markings. fine focus sensitivity 1μm, coarse range 14mm	
Nosepiece	Quintuple, reversed type with positive click stops	
Stage	Low position, mechanical stage, size 180x150mm, moving range 80x50mm.	
	Drop specimen holder(φ110)	
	Specimen holder	
Illumination	Metallurgical reflected illuminator	
	100W halogen lamp house,	
	Electrical parts for 12V 50W halogen lamp	
	12V50W Halogen bulb (Osram)	
	250V/5A Φ5X20 Fuse	
Filter (Φ32)	Power supply 110-240V 50/60Hz	
	Blue filter	
	Green filter	
	Yellow filter	
Polarizing Attachment	Analyzer	
	Polarizer	
Other	Allen screw driver 1.5mm	
	Allen screw driver 2mm	
	Allen screw driver 3mm	
	Metallurgical micrometer(1/100mm)	
	Dust cover	
	Power cord	
Optional Accessories		Item No.
Head	Seidentopf binocular head, 30° inclined, interpupillary distance 52-75mm, diopter adjustable	A53.2702-3
	Seidentopf binocular head, 30° inclined, interpupillary distance 52-75mm, diopter adjustable (For FOV 22mm)	A53.2702-2B
	Seidentopf trinocular head, 30° inclined, interpupillary distance 52-75mm, diopter adjustable (For FOV 22mm)	A53.2704-2A
Eyepiece	WF10x/22mm plan eyepiece, high eyepoint	A51.2701-10C
	WF10x/20mm reticule eyepiece, high eyepoint, diopter adjustable	A51.2703-10A1
Objective	Infinity plan achromatic objective 5x	A52.2706-05A
	Infinity plan achromatic objective 100x	A52.2706-100A
Filter	White Filter (Φ32)	A56.2706-W1
Adapter	1.0x C-Mount	A55.2701-10
	0.5x C-Mount	A55.2701-05
	0.25x C-Mount (suitable to 1/3" CCD)	A55.2701-25
	Digital camera adapter (need to know camera brand)	A55.2702-2
Micrometer	Optical Metallurgical Micrometer	A58.2701





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China