



40X - 400X Metallurgical Polarizing Light Microscope WF10X - 18mm A15.2603

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A15.2603
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Laboratory, Metal
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Magnification: 40X - 400X
- Theory: Polarizing Light Microscope
- Highlight: **student compound microscope ,
compound light microscope**



Product Description

Specifications:

Specification	A15.2603-A	A15.2603-B
Head	Binocular, Inclined 30°, Rotatable 360°, Interpupillary Distance: 50-75mm.	Trinocular, Inclined 30°, Rotatable 360°, Interpupillary Distance: 50-75mm.
Eyepiece	WF10X/18mm	
	WF10X/18mm, with micrometer 0.1mm/Div	
Objective	Infinity Stress Free Plan Objective 4x/0.10	
	Infinity Stress Free Plan Objective 10x/0.25	
	Infinity Stress Free Plan Objective 20x/0.40(S)	
	Infinity Stress Free Plan Objective 40x/0.65(S)	
Nosepiece	Quadruple, Center Adjustable	
Stage	Round stage, Rotatable 360°, Center Adjustable, Division 1°, Vernier Division 6'	
Condenser	Abbe N.A. 1.25, with Iris Diaphragm, Stress Free	
Focusing	Coaxial coarse and fine focusing knob	
Polarizing Set	Analyzer 0-90°, It Can Be Moved Out of the Optical Path for Single Polzrizing Observe	
	built-in Bertrand Llens, Center Adjustable	
	λ Plate, First Class Red	
	λ/4 Plate	
	Quartz Wedge Plate, I-IV Class	
	Poalizer	
Filter	Blue	
Light Source	6V 20W Halogen Lamp, Brightness Adjustable	6V 30W Halogen Lamp, Brightness Adjustable
Optional Accessories		Item No.
Objective	Infinity Stress Free Plan Objective 60x/0.8(S)	A5P.2630-60
Stage	Polarizing Attached Mechanical Stage	A5P.2616
Light Source	Upgrade to 6V 30W	A56.2650
Filter	Green	A56.2616-SG
	Neutral	A56.2616-SN
	Yellow	A56.2616-SA
Adapter	CCD Adapter 0.47x	A55.2601-S47
	CCD Adapter 1.0x	A55.2601-S10





Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China