



## 40x - 1000x Sliding Dual Head Microscope For School A17.1109

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A17.1109
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Application: Lab|educational|school
- Certification: CE|Rohs
- Illumination: Halogen Lamp
- Magnification: 40x - 1000x
- Product Name: Microscope
- Theory: Multi Viewing Microscope
- Highlight: **dual head microscope, Dual View Microscope**



## Product Description

### 40x - 1000x Sliding Dual Head Microscope For School A17.1109

#### Description:

Multi Viewing Microscope, through a multiple optical path, allows 2~6 people observation at the same time. With the features of high quality optical system, effective illumination, LED pointer, the main operator can easily show observation and process to colleagues or students. Coherence, it is widely used in Clinic, Scientific Research and Teaching Demonstration

#### Specifications:

A17.1109 Multi-Viewing Microscope, 2 Position	
Head	Dual Viewing Sliding Binocular Head, 30°Inclined, Interpupillary Distance 55-75mm
Eyepiece	WF10X/18mm, 4 pcs
Objective	Achromatic 4x, 10x, 40x(s), 100x(s, oil)
Nosepiece	Quadruple Revolving Nosepiece
Stage	Double Layer Mechanical Stage 140mm x 140mm, Travel Range 75x50mm
Condenser	Abbe N.A.=1.25 with Iris Diaphragm & Filter
Focusing	Coaxial Coarse & Fine Focusing Dadjustmen Fine Focusing Scale Value 0.002mm
Light Source	Halogen Bulb 6V/20W,AC 220V/110V, Adjustable Brightness, Kohler Illumination



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China