



Gems Diamond Jewelry Microscope Dark Field Microscopes , Rohs CE A24.1205

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A24.1205
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Research, Diamond, Gems
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: LED
- Magnification: 7.5x - 50x
- Theory: Jewelry Microscope
- Highlight: **laboratory binocular microscope,
diamond microscope**



for more products please visit us on cnoec.com

Product Description

Description:

Jewelry Microscope, is special designed microscope on the base of stereo microscope, with jewelry clamp, dark field light source, and more illumination ways, for jewelry and gem inspection and observation.

Specifications:

A24.1205 Jewery Microscope	
Head	Inclined Binocular Head
Eyepiece	WF10x
Objective	0.75x~5x
Transmitted Light	LED Lamp
Working Stage	Jewel Tweezer
	Dark Field Working Stage, With Diaphram

About us :

OPTO-EDU have been focusing in microscope and educational instrument export market for many years. Our customers are from U.S.A, U.K., Russia, Canada, Germany, Denmark, Poland, Sweden, U.A.E. Qatar, Saudia Arabia, Egypt, Argentina, Korea, Thailand, Brazil, and more countries. Concentrating on Optical and Educational products, we are committed to be clients' excellent consultant and supplier for products and services in this field. We strive to be a unique, long-term and excellent supplier for all kinds of microscope from China. We sincerely welcome customers from all over the world contact us at any time!



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China