



24V 6500K Microscope Ring Light Microscope Accessories A56.1211

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.1211
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE, Rohs
- Product Name: Microscope Accessories
- Output Voltage: 24V
- Color Temperature: 6500K
- Brightness Adjustment: 0~100% Adjustable
- Highlight: microscope accessory,
compound microscope parts



More Images



for more products please visit us on cnoec.com

Product Description

Specifications:

Model	A56.1211-24T	A56.1211-48T
LED Bulbs	24 pcs	48 pcs
Output Voltage	24V	
Color Temperature	6500K	
Brightness Adjustment	0~100% Adjustable	
Brightness	At Height of 100mm Around 6500lux	
Working Life	10000 Hours	





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China