



6.5W 460nm Aluminum LED Ring Light , AC 100V - 240V A56.1212

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A56.1212
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE
- Product Name: Microscope Accessories
- Adapter Inner Dia.: $\phi 60\text{mm}$
- Casing Material: Aluminum
- Working Distance: 50-150mm
- Highlight: **compound microscope parts,
parts of the microscope**



More Images



for more products please visit us on cnoec.com

Product Description

Description:

All kinds of Microscope Accessories needed for your microscope, working stage, light source, eyepiece and stage micrometer, digital eyepiece, digital camera, prepare slides, accessories kit, and more.

Specifications:

Model	A56.1212 LED Ring Light, 60 pcs
LED Bulbs	60 pcs
Adapter Inner Dia.	φ60mm
Casing Material	Aluminum
Working Distance	50~150mm
Brightness Control	Bright Controled by Step: 1/4, 2/4. 3/4., 4/4
Brightness	Around 16000 lux, W.D.100mm
Color Temperature	6500K
Working Life	10000 Hours
Wave Length	Blue, 460nm
Power	6.5W
Output Voltage	DC 12V
Output Current	400mA
Input Voltage	AC 100~240V







Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China