



5.0M CMOS LCD Digital microscope camera Microscope Accessories A59.0920

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.0920
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE, Rohs
- Product Name: Microscope Accessories
- Resolution: 5.0M
- Senser: CMOS
- LCD: 3.5" High Clear LCD
- Highlight: **microscope accessory, compound microscope parts**



for more products please visit us on cnoec.com

More Images



Product Description

Description:

A59.0920 3.5: LCD Digital Camera

- ★ 3.5" LCD With Clear Image, Exquisite Outline, Convenient to Carry
- ★ 5.0M Pixels Camera Built-in
- ★ HDMI, AV, USB 2.0 Output All Included for High Resolution Image
- ★ 10M RAM, 4G Kingston SD Card, Large Storage Capacity Achieves Long Working Time
- ★ Can be Mounted on All Kinds of Trinocular Microscope with Standard C-Mount
- ★ Photo/Video Modes For Option, Record Details at Any Time
- ★ Connected To Projector, Image Can Be Magnified For Use of Teaching and Demonstration

Specifications:

A59.0920 LCD Digital Camera Specification	
LCD	3.5" High Clear LCD
Camera	5.0M CMOS Camera Intergrated With SD Card Slot
Sensor	1/2.5"
Resolution	Photo : 2592*1944, 2048*1536, 1024*768 Video : 640*480
White Balance	Automatic White Balance
Image Output	HDMI, AV Interface, USB 2.0
SD Card	4G, Maximum 8G
Working Voltage	External Transformer, 100~240V/5V DC (USB Power Supply)
Support System	Windows 98/ME/2000/XP/VISTA/Win7, Install Relevant Drive Program
Software Funtion	--Date Setting --Language Setting: Chinese, English, Japanese --Manual Adjustment of Exposure Time --Image Previewing (Automatic/Manual) --Auto Power-Off



A59.0920 3.5" LCD Digital Camera



A59.0920 3.5" LCD Digital Camera, Output to HD/Monitor, Computer, Projector, SD Card



Male Tapeworm Strobila

Femal Mosquitoes Mouthparts

Diamond





A59.0920 System Size

**A59.0920 Software has Chinese, English, Japanese Version
Output Include Video, USB2.0 and HDMI, SD Card**



Full Accessories: AV Signal Cable 3.0m, USB 2.0 Cable 2.0m, HDMI Cable 5.0m, Transformer 100~240V/5V DC, SD Card 4G

Connected To Projector, Image Can Be Magnified For Use of Teaching and Demonstration



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China