



1/2 Inch Microscope Digital Cameras Microscope Accessories A59.1003-50B

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.1003-50B
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE
- Product Name: Microscope Accessories
- Resolution: 5.0M
- Sensor: CMOS
- Shutter Type: ERS
- Highlight: **compound microscope parts, parts of the microscope**



More Images



for more products please visit us on cnoec.com

Product Description

Description:

- With high resolution sensor, providing high quality image and widely used in academic and medical field for high precision and resolution image capturing and processing;
- Very convenient to Use, just insert it into Eyepiece Tube or Top Tube of Trinocular Head;
- High sensitivity, low noise, excellent color correction
- USB2.0 live video transfer
- Compatible with Windows 2000/ XP/ Vista/Win 7
- Easy and fast software installation
- User-friendly and sophisticated image processing applications
- Premium optical lens adapter
- Multi-Language software (Chinese, English, German, French, Japanese, Polish, Arabic);
- With measurement function, including angle, length and square.

1/2 Inch Microscope Digital Cameras Microscope Accessories A59.1003-50B

Specifications:

A59.1003-50B Digital Camera		
Image Sensor	1/2 inch 5.0MP CMOS Image Sensor(4:3)	
Active pixels	2592×1944	
Pixel size	2.2μm2.2μm	
Shutter type	Electronic rolling shutter (ERS)	
Frame rate	6fps@2592×1944	
Sensitivity	>0.53 V/lux-sec (550nm)	
Dynamic range	61dB	
SNR	43dB	
Operating temperature	0°C ~ 60°C	
Image Output	USB2.0, Directshow and Twain	
Software	ScopelImage 9.0	
Operating System	Windows 2000/ XP/ Vista/Win7/Win 8/Mac	
A59.1003 Packaging Contents		
Cable	USB2.0 cable	
Software	Software ScopelImage 9.0	
Adapters	Eyepiece Tube Adapter Dia.23.2mm, Ring Φ23.2mm/Φ30.0mm, Ring Φ23.2/Φ30.5mm	
Optional Accessories		
Software	ScopelImage Advanced	Item No.
Stage Micrometer	Calibration slide 0.01mm and 0.1mm	A30.1005







Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China