



CMOS Microscope USB / AV Camera Microscope Accessories A59.1007

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.1007
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE, Rohs
- Product Name: Microscope Accessories
- Sensor: CMOS
- White Balance: Auto
- Connecting Size: Φ 23.2mm
- Highlight: microscope accessory,
compound microscope parts



More Images



for more products please visit us on cnoec.com

Product Description


Description:

- USB & RCA Simulatneous Output
- Excellent Cost Performance

Specifications:

A59.1007 USB/AV Camera		A	B
Digital Camera	1/3" CMOS	●	●
Connecting Size	Φ 23.2mm	●	●
Valid Pixel	720X576	●	●
White Balance	Auto	●	●
Output Mode	USB	●	●
	RCA	●	●
Medium Adapter	0.45X	●	●
Software	ScopelImage 9.0	●	●
System	PAL	●	-
	N	-	●
Ring	Φ 30mm	●	●
Adapter	9V	●	●





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

