



OPTO-EDU A36.5101 LCD USB Stereo Optical Microscope Portable Dual Lens Digital 2.0M 1.3M

Our Product Introduction

for more products please visit us on cnoec.com


Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A36.5101
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




Product Specification

- LCD Screen: "7"" IPS Screen, With Touch Button To Take Photo & Video, One Key Switch Between Main Body Camera & Portable USB Camera"
- Handheld Camera: "1.3M Digital Camra, Zoom Lens 1.5x~25x, Output To 7"" LCD Screen 8 LEDs Light, Brightness Adjustable"
- Main Body Camera: "2.0M Digital Camra, Zoom Lens 5x~35x, TF Card, USB Output Coaxial LED Ring Light, Brightness Adjustable"
- TF Card Output: Support Max 32G, To Take Photo & Video
- USB Output: To Win/Mac, With Professional StrangeView Measuring Software
- Stand: Pole Stand, With Focusing Knob, Size 180x170mm
- Highlight: **usb digital microscope, infinity plan microscope**



Handheld Viewing

More Images



View on Computer

Product Description

7" LCD Stereo + USB Portable Dual Lens Digital Microscope, 2.0M+1.3M

Brand New Design Dual Lens LCD + Portable Digital Microscope For Different Purposes

7.0" LCD, Main Body 2.0M Camera, 5x-35x Zoom For Real Time Clear Watching

Portable USB2.0 1.3M Digital Microscope, 1.5x-25x Zoom For Large Size Objective Watch on LCD

USB Output To Win/Mac, With Professional StrangeView Measuring Software

Support Max 32G TF Card To Take Photo & Video, Power By 5V2A Adapter



A36.5101

7" LCD Stereo + USB Portable Dual Lens Digital Microscope, 2.0M+1.3M



Metal fuselage



Plated film lens



Taking pictures



Connect the computer

LED light source

Touch control

HD

16 million pixels



Adjustable zoom



Super clear quality

Product Details

Metal Zoom Lens 5x~35x

Aviation Aluminum, CNC Carved Molding
High-definition Coating Zoom Lens, The
Focus Is More Precise, The Image Is Clearer



Handheld Viewing

Is More Convenient
Large Size PCB Can Also Be Tested
Quickly



View on Computer

Good Helper For Electronic Repair
Dual-lens Design, More Convenient For
Detection And Storage



LED Illumination

Stepless Dimming Of The Main
Lens Hold The Lens And Touch The
Button To Dim

Metal body

Aviation Aluminum, Sturdy And
Durable



Switch &
Handheld Camera
Or Main Body
Camera Switch

Menu



MENU



OK



Take Photo
& Menu Up

Take Video &
Menu Down

Photo Album &
Menu Confirmation

TF Card


USB

Power Supply

Work Indicator




Practical Application



Clock & Electronic Precision Repair





Lapis Lazuli

Emerald Rough Stone

Gem Detection




Detection Of Various Currencies & Commemorative Coins & Banknotes & Stamps



Textile Testing

PC Measuring Software





A36.5101 7" LCD Stereo + USB Portable Dual Lens Digital Microscope, 2.0M+1.3M	
LCD Screen	7" IPS Screen, With Touch Button To Take Photo & Video, One Key Switch Between Main Body Camera & Portable USB Camera
Handheld Camera	1.3M Digital Camra, Zoom Lens 1.5x~25x, Output to 7" LCD Screen 8 LEDs Light, Brightness Adjustable
Main Body Camera	2.0M Digital Camra, Zoom Lens 5x~35x, TF Card, USB Output Coaxial LED Ring Light, Brightness Adjustable
TF Card Output	Support Max 32G, To Take Photo & Video
USB Output	To Win/Mac, With Professional StrangeView Measuring Software
Stand	Pole Stand, With Focusing Knob, Size 180x170mm
Size	180*170*322mm
Power Supply	DC 5V2A Adapter, Can Be Powered By USB Power Bank For Outdoor Use
Package	Color Box + Foam, 1 pc/1 Box, With Software Disc, DC Adapter, Instruction Manual



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



sale@optoedu.com



cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China