



OPTO-EDU A59.4231 2K 38M Portable Hd Microscope Camera

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.4231
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Sensor: Panasonic CMOS 21M, 1/2.3"
- Pixel Size: 1.335 μ m X 1.335 μ m
- Output: HD(Type A), USB(Type B),TF Card
- Photo Resolution: 38M
- HD Support: 4:3, 16:9, 16:10
- White Balance: Auto/Manual, RGB Adjustable Separately
- Highlight: **digital trinocular microscope, usb digital microscope**



More Images




Product Description


ccd camera for microscope

OPTO-EDU A59.231 38mp Digital Eyepiece Camera For Microscope With Electronic Usb HD USB Digital Camera For Trinocular Microscope

Our Strength

25+ Years Professional Experience Know Microscopes Best
150+ Microscope & Accessory Manufacturers From China
200+ Hot Sale & Newest Update Microscope Model Every Month
720+ Customers From All Over World Keep Rising
1200+ Microscope Models To Create One-Stop Purchase Platform
OPTO-EDU As Famous Brand For Microscope In China & USA
Top Seller For Microscope on Alibaba Recommend To 99% User
Most Professional Sale Team Supply Best Microscope Solution







A59.4231 2K HDMI Digital Camera, 38.0M	
Sensor	Panasonic CMOS 21M, 1/2.3"
Pixel Size	1.335 μ m \times 1.335 μ m
Output	HDMI(Type A), USB(Type B),TF Card
Photo Resolution	38M
Video Resolution	2K@30fps, 1080p@60fps USB2.0 Output (No Driver Need) 1920*1080@30 1280*720@30 640*480@30 HDMI 1920*1080 60fps, Output With USB at same time
HDMI Support	4:3, 16:9, 16:10
TF Card	Support Max 64GB
White Balance	Auto/Manual, RGB Adjustable Separtely
Exposure	Auto/Manual, Adjustable Exporsure
Cross Line	8 Lines, Color, Width, Position Adjustable; 2 Lines With Scale
Menu	Sharpness, Contrast, Saturation Adjustable
Black/White, Negative Film	Support
Mirror Image	Support
Language	Chinese/English
Digital Zoom Support	With Remote Control Photo, Video
Input	DC 12V (Support 5-19V), 55-21 Socket
Power	Consumption Max 0.15W (12V,0.12A)
Adapter	C/CS Mount
Net Weight	180g
Gross Weight	400g
Size	74*74*35mm
Package Size	160*100*80mm



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China