



## Research USB Handheld Digital Microscope Digital Camera Microscopes A34.5001

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.5001
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Model Number: A34.5001
- Application: Research, Learning, Home
- Certification: CE
- Magnification: 20x-300x
- Pixel: 5 Mega
- Theory: Handheld Digital Microscope
- Highlight: **Portable Digital Microscope,  
Hand held Digital Microscope**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Research USB Handheld Digital Microscope Digital Camera Microscopes A34.5001

#### Description:

- USB Digital Microscope with professional Stand
- 5M
- 20x-300x Magnification
- Measurement
- Mac compatible

#### Specifications:

A34.5001 USB Digital Microscope	
Image Sensor:	5 Mega Pixels (true resolution)
Still capture resolution	5M, 3M, 1.3M, VGA, QVGA
Video capture resolution	5M, 3M, 1.3M, VGA, QVGA
Focus Range	Manual focus from 10mm to 50mm
Frame Rate	Max 30f/s under 600 Lus Brightness
Magnification Ratio	20x to200x and 300x
Video format	AVI
Photo format	JPEG or BMP
Light source	8 LED (adjustable by control wheel)
PC interface:	USB2.0
Operation system:	Windows XP/Vista/Win7 Mac OS X 10.6 or above
OSD language	English, German, Spanish, French, Russian
Bundle software	MicroCapture with measurement
Product Size	110mm (L) x 33mm (R)
Weight	90g
Warranty	1 year









**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China