



## Silver 1.3M CMOS USB Handheld Digital Microscope For Learning A34.5501

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A34.5501
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Model Number: A34.5501
- Application: Learning, Home
- Certification: CE
- Magnification: 20x-350x
- Pixel: 1.3M CMOS
- Theory: Handheld Digital Microscope
- Highlight: **Portable Digital Microscope,  
Hand held Digital Microscope**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Silver 1.3M CMOS USB Handheld Digital Microscope For Learning A34.5501

#### Description:

Digital Hand-held Microscope, designed in mini size, portable outlook, is the best tool for you to learn and find micro world around you. You can use it to see almost everything in micro structure on a computer, or even your iPad, iPhone, Android mobile devices.

#### Specifications:

A34.5501 USB Digital Microscope	
Magnification	20x-350x
Pixel	1.3M CMOS
Resolution	640*480 (can adjust from 640*480 to 1280*1024)
Video format	AVI
Photo format	JPEG or BMP
Light source	Top:4pcs LED (adjustable brightness) Bottom:1 Led (adjustable brightness)
PC interface:	USB
	Manual instruction
Packing list	CD(Microscope operation, image edit software and driver software) Microscope tools: Forceps, probes, slides, samples, etc.





**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ [sale@optoedu.com](mailto:sale@optoedu.com)

🌐 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China