



## OPTO-EDU A59.3520 2048x1536 All In One 12.0m Opto Edu Microscope Lcd Screen

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.3520
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Theory: Digital Camera
- Resolution: 2048x1536
- CPU: 1.8GHz Four Core RK3288 CPU
- WIFI: Support WIFI & Bluetooth
- Product Size: Pad 18.2\*23.7\*1.1cm
- Operation System: Android 5.1 System, Multi Language
- Highlight: **digital trinocular microscope, usb digital microscope**



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### OPTO-EDU A59.3520 9.7" LCD Android Pad Digital Camera For Microscope All-in-One, 12.0M

5.0M, 9.7" LCD, Resolution 2048\*1536, Android 5.1, CPU 4 Core 1.8G, RAM 2G ROM 8G, HDMI+TF Card Output, Touch Screen USB Mouse Control, Installed Measure Software, Chinese/English Manual, Size 18.2\*23.7\*1.1cm, Power 100~240V, 50/60Hz

9.7" LCD, Resolution 2048\*1536

5.0M 1/2.5" CMOS Digital Camera Built-in

Android 5.1, CPU 4 Core 1.8G, RAM 2G ROM 8G

HDMI+TF Card Output, Touch Screen USB Mouse Control

Installed Measure Software, Chinese/English Manual



### 9.7" Android Pad + 5.0M Camera

Sharp Image | All-in-One | Photo & Video



for more products please visit us on [cnoec.com](http://cnoec.com)

Our Product Introduction

## Product Details




HDMI And TF-Card Interface

Mini USB2.0 & USB2.0 Interface, It Has Mike And Headphone Function

## Measuring Software Included

Touch Screen | Full Range Measure | 90 Languages



## Product Applications



A59.3520 9.7" LCD Pad Digital Camera	
LCD	9.7" HD IPS Touch Screen With Backlight
Camera	5.0M 1/2.5" CMOS Digital Camera
Resolution	2048x1536 @15FPS
Pixel	2.2umx2.2um
Data Storage	TF Card or Built-in Memory
Data Output	USB2.0, HDMI, TF Card
Photo	5.0M .JPEG
Video	1080p Video Record
System	Android 5.1 System
	90 Languages
	Support Mouse + Keyboard Control
	Support WIFI
Function	Support Bluetooth
	S-Eye Software & Instruction Manual Installed
	Mark: Point, Crosshair, Coordinate, Text
	Length Measure: Straight Line Length, Broken Line Length, Curve Length, Parallel Line Distance, Point To Line Distance
	Distance Between 2 Circle Centers: Radius Fixed Circle, Two-point Fixed Circle, Three-point Fixed Circle
	Geometric Measure: Line Segment Length, Radius Fixed Circle, Two-point Fixed Circle, Three-point Fixed Circle, Concentric Circle
	Geometric Area Measure: Polygon, Square
Power	DC 12V/2A Power Adaptor, AC 100~240V
Size	Pad 238x182x12mm, Camera 40x39mm
Package	Packing Box 10x24x29cm, 0.6Kg, 1 pc/Box



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China