



400X / 100X Laboratory Compound Optical Microscope , Phase Contrast Microscopes A12.1007

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1007
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory, Hobby, Home
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp, LED
- Max Magnification: 1000X
- Product Name: Compound Optical Microscope, Microscope
- Highlight: **binocular compound microscope , optical compound microscope**



More Images



for more products please visit us on cnoec.com

Product Description

Laboratory 1000X Compound Optical Microscope , Rohs CE A12.1007

Specifications:

A12.1007 Biological Microscope		
Head	Siedentopf binocular head, inclined at 30°,360° Rotatable	
Eyepiece	Wide field eyepiece WF10x/ 18	
Objective	Achromatic Objective 4x	
	Achromatic Objective 10x	
	Achromatic Objective 40x	
	Achromatic Objective 100x	
Nosepiece	Backward quadruple nosepiece	
Stage	Double Layers Mechanical Stage 140×140mm/ 75×50mm	
Focusing	Coaxial coarse & fine adjustment system, range 28mm, fine division 0.002mm	
Condenser	Abbe condenser N.A.1.25 with Iris diaphragm & filter,rack & pinion adjustment	
Filter	Green Filter	
Illumination		
Optional Accessories		Item No.
Head	Monocular head, inclined at 30°,360° Rotatable	A53.1002-3
	Seidentopf trinocular head, inclined at 30°,360° Rotatable	A53.1004-2
Eyepiece	Wide Field Eyepiece WF16x/ 11	A51.1002-16
Objective	Achromatic objective 20x	A52.1001-20
	Achromatic objective 60x	A52.1001-60
	Plan Achromatic Objective 4x	A52.1006-4
	Plan Achromatic Objective 10x	A52.1006-10
	Plan Achromatic Objective 20x	A52.1006-20
	Plan Achromatic Objective 40x	A52.1006-40
	Plan Achromatic Objective 60x	A52.1006-60
	Plan Achromatic Objective 100x	A52.1006-100
Nosepiece	Backward quintuple nosepiece	A54.1001-5M
Illumination	Halogen lamp 6V/ 20W, brightness adjustable	A56.1003-6V20
	Kohler illumination	A56.1005-3
Video	1X CMount	A55.1002-10
Attachment	0.5X CMount	A55.1002-05
Epi-fluorescent Attachment	XY-2 EPI-fluorescent Attachment	A5F.1001
Attachment	FL-LED EPI-fluorescent Attachment	
Phase	Simple phase contrast set	A5C.1001
Contrast	Sliding phase contrast set	A5C.1002-C3
Attachment	Turret phase contrast set	A5C.1003-B1
Polarization Set	Simple polarization set	A5P.1001
Dark Field	Dark field condenser (Dry) N.A. 0.9	A5D.1010-BD
Condensor	Dark field condenser (Oil) N.A. 1.25-1.36	A5D.1010-BI







Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China