



Student Compound Optical Microscope

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1018
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School, Laboratory
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp, LED
- Max Magnification: 1000X
- Product Name: Compound Optical Microscope
- Highlight: **binocular compound microscope, wide field microscope**



for more products please visit us on cnoec.com

More Images



Product Description

Student Compound Optical Microscope Halogen Lamp Microscopes A12.1018

Specifications:

A12.1018 Biological Microscope		
Head	Siedentopf trinocular viewing head	
Eyeiece	Extra wide field eyepiece EW10x/ 22	
Objective	Infinite plan achromatic objective 4x	
	Infinite plan achromatic objective 10x	
	Infinite plan achromatic objective 40x	
	Infinite plan achromatic objective 100x	
Nosepiece	Backward quadruple nosepiece	
Focusing	Coaxial coarse & fine adjustment, range 24mm, fine division 0.001mm,coarse stroke 37.7mm	
Stage	Double Layers Mechanical Stage 185×142/ 75×55mm	
Condenser	Swing-out Condenser NA 0.9/ 0.25	
Filter	Blue filter	
	Green filter	
Illumination	External Illumination, Aspherical Collector, Halogen Lamp 6V/ 30W	
Optional Accessories		
	Item No.	
Head	Seidentopf Binocular Head	A53.1003-7
Objective	Infinity plan achromatic objective 20x	A52.1006-20
	Objectives for Dark Field 100x	A5D.1002-3
Nosepiece	Backward sextuple nosepiece	A54.1001-6
Illumination	LED	A56.1004-3
	24V/100W Power Supply	A56.1011-800
Photo Attachment		A55.1001
Video Adapter	C Mount 1x	A55.1002-10
	C Mount 0.5x	A55.1002-05
Dark Field Attachment	Dark field condenser (Dry) N.A. 0.9	A5D.1010-BD
	Dark field condenser (Oil) N.A. 1.25-1.36	A5D.1010-BI
Phase Contrast Attachment	Turret phase contrast set	A5C.1003-A1
Polarization Set		A5P.1003
Epi-fluorescent Attachment	FL-800 EPI-fluorescent Attachment	A5F.1002
	FL-LED EPI-fluorescent Attachment	





www.cnoec.com



www.cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China