



## 20X Stereo Optical Microscope

Our Product Introduction

for more products please visit us on cnoec.com

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1204
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Magnification: 20X
- Applications: Medical, Clinic
- Certification: CE, Rohs
- Drawtube: Binocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: stereo dissecting microscope, stereo binocular microscopes



### More Images



### Product Description

#### 20X Clinic Stereo Optical Microscope Low Magnification Microscopes A22.1204

##### Description:

Stereo Optical Microscope, designed for low magnification observation of a sample, generate three-dimensional, laterally precise images, to see where a three-dimensional view is needed for analyzing the detail. It is ideal for observing large biological samples or analysis of rough materials surfaces, widely used in industry production line, coin & jewelry inspection and etc.

##### Specifications:

A22.1204 Specification				
Model	Eyeiece	Objective	Working Stage	Illuminator
1A	WF10X	Inserted 2x	White/Black	No
1B	WF10X	Inserted 2x	White/Black	Incident
1C	WF10X	Inserted 2x	Frosted Glass	Incident & Transmitted
Optical Accessories				Item No.
Eyepiece	WF5X			A51.1201-05
	WF10X			A51.1201-10
	WF15X			A51.1201-15
	WF20X			A51.1201-20
	1x			A52.1280-01

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

Objective	2x					A52.1280-02			
	3x					A52.1280-03			
	4x					A52.1280-04			
	6x					A52.1280-06			
Technical Data									
Eyepiece									
Objec	5X		10X		15X		20X		Working Distance (mm)
	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Magnification	Field of View (mm)	
1X	5X	20.0	10X	20.0	15X	15.0	20X	10.0	57
2X	10X	10.0	20X	10.0	30X	7.5	40X	6.5	80
3X	15X	7.0	30X	6.7	45X	5.0	60X	4.3	61
4X	20X	5.0	40X	5.0	60X	3.7	80X	3.2	57
6X	30X	4.0	60X	3.3	90X	2.5	120X	2.1	66





Opto-Edu (Beijing) Co., Ltd.

0086 13911110627 sale@optoedu.com cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China