



## Medical Stereo Optical Microscope

Our Product Introduction

for more products please visit us on cnoec.com

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A22.1205
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Magnification: 20X
- Applications: Gem, Jewelry, Medical
- Certification: CE, Rohs
- Drawtube: Binocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: stereo dissecting microscope, stereo binocular microscopes



### More Images



### Product Description

#### Medical Stereoscopic Microscope Cordless Microscopes , Rohs CE A22.1205

#### Description:

Stereo Optical Microscope, designed for low magnification observation of a sample, generate three-dimensional, laterally precise images, to see where a three-dimensional view is needed for analyzing the detail. It is ideal for observing large biological samples or analysis of rough materials surfaces, widely used in industry production line, coin & jewelry inspection and etc.

#### Specifications:

A22.1205 Specification				
Model	Eyeiece	Objective	Plate	Illminator
2A	WF10X	Inserted 2x	White/Black	No
2B	WF10X	Inserted 2x	White/Black	Incident
2C	WF10X	Inserted 2x	Frosted Glass	Incident & Transmitted
Optical Accesies				Item No.
Eyepiece	WF5X			A51.1201-05
	WF10X			A51.1201-10
	WF15X			A51.1201-15
	WF20X			A51.1201-20

Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)

Objective	1x							A52.1280-01	
	2x							A52.1280-02	
	3x							A52.1280-03	
	4x							A52.1280-04	
	6x							A52.1280-06	
Technical Data									
Eyepiece									
	5X		10X		15X		20X		
Objective	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Magnification	Field of View (mm)	Working Distance (mm)
1X	5X	20.0	10X	20.0	15X	15.0	20X	10.0	57
2X	10X	10.0	20X	10.0	30X	7.5	40X	6.5	80
3X	15X	7.0	30X	6.7	45X	5.0	60X	4.3	61
4X	20X	5.0	40X	5.0	60X	3.7	80X	3.2	57
6X	30X	4.0	60X	3.3	90X	2.5	120X	2.1	66





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627    ✉ sale@optoedu.com    🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China