



Portable Biological Microscope

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1123
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School, Hobby, Home
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: LED
- Max Magnification: 1000X
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **binocular biological microscope, biological compound microscope**



More Images



Product Description

Portable Cordless Biological Microscope 1000X For School / Home A11.1123

Description:

We have various Student Biological Microscope for your choice for student educational purpose. Designed with high quality optical system, steady and heavy materials, we can supply a real optical microscope for all students at very good prices. Sold to world market in large quantity, these student microscopes has become most popular tools in primary, middle and high school.

Specifications:

A 11.1123 Economical Binocular Biological microscope

This scope is an ideal school microscope, it is with LED lighting and fully portable with no need for messy classroom wires! Arm with handle hole is easy to carry

Our Product Introduction

for more products please visit us on cnoec.com

A11.1123 Economical Binocular Biological Microscope Specification	
Head	Sliding Binocular Head, 45°Inclined, 360°Rotatable
Eyeiece	WF10X with pointer Achromatic 4X
Objective	Achromatic 10X Achromatic 40X, Spring Achromatic 100X, Spring/Oil
Nosepiece	Quadruple Revolving Nosepiece
Stage	Double layer Mechanical Stage ,Size 100x90mm
Condenser	NA1.25 Condenser with Iris Diaphragm
Focusing	Coaxial Coarse and Fine Focusing Knobs
Light Source	Top and Bottom LED light with rechargeable function (if shipped by air, battery will not be included) Power supply adapter




www.cnoec.com



www.cnoec.com



www.cnoec.com





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China