



Fluorescent Jewelry Microscope

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A24.1202
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Research
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: Halogen Lamp
- Magnification: 7X - 40X
- Theory: Jewelry Microscope
- Highlight: **laboratory binocular microscope,
gemological microscope**



Product Description

Fluorescent Jewelry Microscope With 360 Rotatable Head , CE A24.1202

Specifications:

A24.1202 Jewery Microscope			
Head	Inclined Binocular Head,360° Rotatable		
Eyepiece	SWF10x		
Objective	1x~4x	•	
	0.7x~3.6x		•
Incident Light	Halogen Lamp 12V/10W, Brightness Adjustable		
Transmitted Light	Fluorescent 10W, Brightness Adjustable		
Working Stage	Jewel Tweezer		
	Dark Field Working Stage, With Diaphram		
A24.1202 Jewery Microscope , Optional Accessories			
Eyepiece	WF5x		E51.1221-05
	WF10x		E51.1221-10
	WF15x		E51.1221-15
	WF20x		E51.1221-20



Our Product Introduction

for more products please visit us on cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China