



40x - 1000x Student Biological Microscope Fine Adjustment Knob Microscope A11.1010

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1010
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student|school
- Certification: CE|Rohs
- Drawtube: Binocular|Monocular
- Illumination: Halogen Lamp
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
binocular biological microscope**

for more products please visit us on cnoec.com

Product Description

40x - 1000x Student Biological Microscope Fine Adjustment Knob Microscope A11.1010

- 195mm Achromatic Objectives, Mechanical Tube Length 160mm
- WF10x/18mm, WF16x/11mm Wide Field Eyepiece
- Monocular, Binocular, Trinocular Head Optional
- Low Position Coaxial Coarse & Fine Focus System
- Double Layer Mechanical Stage

| A11.1010 Biological Microscope | | A1 | A2 | A3 | B |
|--|--|----|----|----|-------------|
| Magnification | 40x~1000x | | | | |
| Head | Monocular Head, 45° Inclined, 360° Rotatable | o | o | o | |
| | Seidentopf Binocular Head, 30° Inclined, 360° Rotatable, Interpupillary Distance 55~75mm | | | | o |
| Eyepiece | Wide Field WF10x/18mm | | | | |
| Nosepiece | Backward Triple Nosepiece | o | o | | |
| | Backward Quadruple Nosepiece | | | o | o |
| Objective | Achromatic 4x/0.10, 10x/0.25, 40xs/0.65 | o | o | | |
| | Achromatic 4x/0.10, 10x/0.25, 40xs/0.65, 100xs/1.25(Oil) | | | o | o |
| Stage | Plain Stage, Moving Range 60*29mm | o | o | | |
| | Mechanical Stage, Moving Range 60*30mm | | | o | o |
| Focusing | Coaxial Coarse & Fine Focusing | | | | |
| Condenser | Disc Diaphragm | o | | | |
| | N.A.1.25 Abbe Condenser, Iris Diaphragm | | o | o | o |
| Light Source | Plan-Concave Mirror | o | | | |
| | Built-in 6V20W Halogen Lamp, Brightness Adjustable | | o | o | o |
| A11.1010 Biological Microscope Optional Accessories | | | | | |
| Eyepiece | P16x | | | | A51.1002-16 |



A11.1010-B Binocular Biological Microscope

A11.1010-A1,A2,A3



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China