



Coarse Fine Adjustment Microscope Trinocular For Student A11.1008

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1008
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, Laboratory, Hobby
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: Halogen Lamp
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
binocular biological microscope**



Product Description

Coarse Fine Adjustment Microscope Trinocular For Student A11.1008

Specifications:

A11.1008 Biological Microscope		A	B
Head	Sliding Binocular head	●	●
	Anti-mould		●
Eyepiece	Wide Field Eyepiece WF10x	●	●
Objective	Achromatic objective 4x	●	●
	Achromatic objective 10x	●	●
	Achromatic objective 40x	●	●
	Achromatic objective 100x	●	●
Nosepiece	Quadruple nosepiece	●	●
Stage	Double Layer Mechanical Stage 140x140mm/ 75x50mm	●	●
Focusing	Coaxial Coarse & Fine knob, Fine division 0.002mm, Coarse stroke 37.7per rotation	●	●
	Rubber cover for focus adjustment knob		●
Condenser	Abbe N.A. 1.25 with Iris Diaphragm & Filter	●	●
Illumination	Built-in Illumination, Halogen Lamp 100v-240W /6V-20W	●	●
	Plan-concave Mirror	●	●
Optional Accessories		Item No.	
Head	Sliding trinocular head	A53.1003-3	
	Seidentopf binocular head	A53.1004-1	
Eyepiece	Wide Field Eyepiece P16x	A51.1002-16	
Nosepiece	Quintuple nosepiece	A54.1001-5M	



A12.1008-A Student Biological Microscope



A12.1008-B Student Biological Microscope



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China