



400x Monocular Biological Microscope With Electric Light A11.1101

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1101
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School
- Certification: CE
- Drawtube: Monocular
- Max Magnification: 400x
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **binocular biological microscope,
biological compound microscope**



More Images



Product Description

400x Monocular Biological Microscope With Electric Light A11.1101

Specifications:

A11.1101 Biological microscope		-A	-B	-C
Head	Monocular Head, Straight Tube	●	●	●
Eyepiece	WF10X with the Pointer	●	●	●
Objective	Achromatic 4X	●	●	●
	Achromatic 10X	●	●	●
	Achromatic 40X, Spring	●	●	●
Nosepiece	Triple Revolving Nosepiece	●	●	●
Stage	Square Stage with Specimen Clips, Size 84x87mm	●	●	●
Condenser	Single Lens with Disc Diaphragm	●	●	●
Focusing	Coarse Focusing Knobs	●	-	-
	Separate Coarse and Fine Focusing Knobs	-	●	●
Light Source	Plan-concave Reflective Mirror	●	●	-
	Electrical Light Source 220V/15W or 110V/15W	-	-	●
Optional Accessories				
Eyepiece	WF16X			

Description:

We have various Student Biological Microscope for your choice for student educational purpose. Designed with high quality optical system, steady and heavy materials, we can supply a real optical microscope for all students at very good prices. Sold to world market in large quantity, these student microscopes has become most popular tools in primary, middle and high school.

About us :

OPTO-EDU have been focusing in microscope and educational instrument export market for many years. Our customers are from U.S.A, U.K., Russia, Canada, Germany, Denmark, Poland, Sweden, U.A.E. Qatar, Saudia Arabia, Egypt, Argentina, Korea, Thailand, Brazil, and more countries. Concentrating on Optical and Educational products, we are committed to be clients' excellent consultant and supplier for products and services in this field. We strive to be a unique, long-term and excellent supplier for all kinds of microscope from China. We sincerely welcome customers from all over the world contact us at any time!





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China