



Achromatic Wide Field Microscopes With Electric Light Source CE A11.1102

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1102
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: School, Laboratory
- Certification: CE, Rohs
- Drawtube: Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **binocular biological microscope,
biological compound microscope**



More Images



Product Description

Achromatic Wide Field Microscopes With Electric Light Source CE A11.1102

Specifications:

A11.1102 Biological microscope		-A	-B	-C
Head	Monocular Head,45° inclined,360° Rotatable	●	●	●
Eyepiece	WF10X with the pointer	●	●	●
Objective	Achromatic 4X	●	●	●
	Achromatic 10X	●	●	●
	Achromatic 40X, Spring	●	●	●
Nosepiece	Triple Revolving Nosepiece	●	●	●
Stage	Square Stage with Specimen Clips, Size 84x87mm	●	●	●
Condenser	Single Lens with Disc Diaphragm	●	●	●
Focusing	Coarse Focusing Knobs	●	-	-
	Separate Coarse and Fine Focusing Knobs	-	●	●
Light Source	Plan-concave Reflective Mirror	●	●	-
	Electrical Light Source 220V/15Wor110V/15W	-	-	●
Optional Accessories				
Eyepiece	WF16X			
	20X			
Objective	60X(S)			
	100X(S)Oil			
Light Source	LED			



www.cnoec.com



A11.1102-B

A11.1102-C



Opto-Edu (Beijing) Co., Ltd.

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China