



Disc Diaphragm Biological Microscope Monocular , 40X - 400X A11.1105

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1105
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student
- Certification: CE
- Drawtube: Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **binocular biological microscope,
biological compound microscope**



Product Description

Description:

We have various Student Biological Microscope for your choice for student educational purpose. Designed with high quality optical system, steady and heavy materials, we can supply a real optical microscope for all students at very good prices. Sold to world market in large quantity, these student microscopes has become most popular tools in primary, middle and high school.

Specifications:

A11.1105 Student Microscope	
Head	Monocular,Inclined at 45°,360° Rotatable
Eyepiece	WF10X
Objective	Achromatic 4x
	Achromatic10x
	Achromatic 40x(s)
Stage	Plain Stage with Slides Clips, 110mmx105mm
Condenser	Single Lens with Disc Diaphragm
Focusing	Separated Coarse & Fine Focusing
Light Source	Incandescent Lamp,220V/20Wor 110V/20W
Optional Accessories	
Eyepiece	WF16X
	P16X
	20X
Objective	60X(S)
	100(S/O)
Light Source	LED

About us :

OPTO-EDU have been focusing in microscope and educational instrument export market for many years. Our customers are from U.S.A, U.K., Russia, Canada, Germany, Denmark, Poland, Sweden, U.A.E. Qatar, Saudia Arabia, Egypt, Argentina, Korea, Thailand, Brazil, and more countries. Concentrating on Optical and Educational products, we are committed to be clients' excellent consultant and supplier for products and services in this field. We strive to be a unique, long-term and excellent supplier for all kinds of microscope from China. We sincerely welcome customers from all over the world contact us at any time!



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China