



Sliding Binocular Biological Microscope For Hobby A11.1106

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1106
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Hobby, Home
- Certification: CE, Rohs
- Drawtube: Binocular, Monocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
biological compound microscope**



More Images



Product Description

Specifications:

A11.1106 Biological Microscope		-A	-B	-C	-D
Head	Monocular Head,45° inclined,360° Rotatable	●	●	●	-
	Sliding Binocular Head,45° inclined,55mm-75mm	-	-	-	●
Eyepiece	WF10X/18mm	●	●	●	●
Objective	Achromatic 4X	●	●	●	●
	Achromatic 10X	●	●	●	●
	Achromatic 40X, Spring	●	●	●	●
	Achromatic 100X, Spring	-	-	●	●
Nosepiece	Triple Revolving Nosepiece	●	●	-	-
	Quadruple Revolving Nosepiece	-	-	●	●
Stage	Plain Stage with Specimen Clips, Size120x120mm	●	●	-	-
	Mechanical Stage,60mmx30mm	-	-	●	●
Condenser	N.A.=0.65, with Disc Diaphragm	●	-	-	-
	N.A.=0.65, with Iris Diaphragm	-	●	-	-
	Abbe Condenser N.A.=1.25, with Iris Diaphragm & Filter	-	-	●	●
Focusing	Seperate Coarse & Fine Focusing Knobs	●	●	●	●
Light Source	Incandescent Lamp220V/20W or 110V/20W or Fluorescent	●	●	-	-
	Halogen Bulb 6V/20W,AC 220V/110V	-	-	●	●
Optional Accessories					
Eyepiece	WF16X				
	P16X				
Objective	Achromatic 20X				
	Achromatic 60X(S)				
Light Source	LED				







Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China