



## 1600x Separate Coarse Biological Microscope For School A11.1107

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1107
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Student, School
- Certification: CE, Rohs
- Drawtube: Monocular
- Max Magnification: 1600x
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,  
binocular biological microscope**



### More Images



for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Specifications:

A11.1107 Biological Microscope		-A	-B	-C
Head	Monocular Head, Straight Tube	●	●	●
Huygens Eyepiece	H5	-	●	-
	H6	-	-	●
	H10	●	●	●
	H15	-	-	●
	H16	●	●	-
Objective	Achromatic 4X	●	-	-
	Achromatic 5X	-	-	●
	Achromatic 10X	●	●	●
	Achromatic 40X, Spring	●	●	-
	Achromatic 45X	-	-	●
	Achromatic 100X, Spring, Oil	-	●	-
Nosepiece	Triple Revolving Nosepiece	●	●	●
Stage	Plain Stage with Disc Diaphragm	●	-	-
	Moving Stage with Scale, Moving Range: 60mmx30mm	-	●	-
	Round Stage with Iris Diaphragm,	-	-	●
Focusing	Separate Coarse and Fine Focusing Knobs	●	●	●
Light Source	Reflective Mirror	●	●	●
Optional Accessories				
Eyepiece	H12.5X			
Light Source	Electric Light Source 220V/15W or 110V/15W			



www.cnoec.com



www.cnoec.com



www.cnoec.com



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China