



School Biological Microscope Kohler Illumination Microscopes A11.1109

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1109
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: School, Laboratory
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Max Magnification: 1000X
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
binocular biological microscope**



for more products please visit us on cnoec.com

Product Description

School Biological Microscope Kohler Illumination Microscopes A11.1109

Specifications:

A11.1109 Biological Microscope		-A	-B
Head	Compensation Free Binocular,30°inclined,55mm-75mm	●	-
	Compensation Free Trinocular,30°inclined,55mm-75mm	-	●
Eyepiece	WF10/18mm	●	●
Objective	Achromatic 4X	●	●
	Achromatic 10X	●	●
	Achromatic 40X, Spring	●	●
	Achromatic 100X, Spring	●	●
Nosepiece	Quadruple Revolving Nosepiece	●	●
Stage	Double Layer Mechanical,Size 140mmX140mm,Moving Range:70mmx50mm	●	●
Condenser	ABBE N.A.1.25Condenser with Iris Diaphragm & Filter	●	●
Focusing	Coaxial Coarse and Fine Focusing Scale Value 0.002mm	●	●
Illumination	Kohler Illumination	●	●
Light Source	Halogen Bulb 6V/20W,AC 220V/110V, Adjustable Brightness	●	●

Optional Accessories	
Eyepiece	WF16X
	WF20X
	P16X
	Achromatic 20X
	Achromatic 60X, Spring
	Semi-plan Objective 4X
	Semi-plan Objective 10X
Objective	Semi-plan Objective 20X
	Semi-plan Objective 40X(S)
	Semi-plan Objective 60X(S)
	Semi-plan Objective 100X(S)
	Plan Objective 4X
	Plan Objective 10X
	Plan Objective 20X
	Plan Objective 40X(S)
	Plan Objective 60X(S)
	Plan Objective 100X(S)
Dark Field	Darkfield
Photo Adapter	Photo Adapter
CCD Adapter	0.5x
	0.57x
	0.75x



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China