



LED Disc / Iris Diaphragm Microscope With 360°Rotatable Head A11.1112

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1112
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School
- Certification: CE, Rohs
- Drawtube: Binocular, Monocular, Trinocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
binocular biological microscope**



More Images



for more products please visit us on cnoec.com

Product Description

LED Disc / Iris Diaphragm Microscope With 360°Rotatable Head A11.1112

Description:

This microscope is suitable for school student use

- Full Coated Optical System.
- ccurate & Smooth Focusing System with Built-in Slip Clutch and Tension Adjustment
- safty Rack Stop Prevents Damage to Slides & Objectives
- Disc or Iris Diaphragms Optional
- Well Painted Metal Base & Stand

Specifications:

A11.1112 Biological Microscope		-A	-B	-C
Head	Monocular Head,45° inclined,360° Rotatable	●	-	-
	Compensation Free Binocular Head,30° inclined,55mm-75mm	-	●	-
	Compensation Free Trinocular Head,30° inclined,55mm-75mm	-	-	●
Eyepiece	WF10X/18mm	●	●	●
Objective	Achromatic 4X	●	●	●
	Achromatic 10X	●	●	●
	Achromatic 40X, Spring	●	●	●
	Achromatic 100X, Spring,Oil	●	●	●
Nosepiece	Quadruple Revolving Nosepiece	●	●	●
Stage	Mechanical Stage with Scale,Size:120mmx125mm Moving Range:60mmx30mm	●	●	●
Focusing	Seperate Coarse & Fine Focusing Knobs	●	●	●
Light Source	Halogen Bulb 6V/20W,Brightness Adjustble	●	●	●
Optional Accessories				
Eyepiece	WF16X			
	WF20X			
	P16X			
Objective	Achromatic 20X			
	Achromatic 60X(S)			
	Semi-plan Achromatic Objective			
	Plan Achromatic Objective			
Light Source	Reflect Mirror			
	Fluorescent Lamp 220V/5W or 110V/5W			
	LED			
CCD Adapter	CCD Adapter			



www.cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China