



Monocular / Trinocular Biological Microscope 400X with Halogen Bulb A11.1116

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1116
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: School, Laboratory
- Certification: CE, Rohs
- Drawtube: Binocular, Monocular, Trinocular
- Illumination: LED
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
biological compound microscope**



More Images



Product Description

Specifications:

A11.1116 Biological Microscope		-A	-B	-C
Head	Monocular Head,45° inclined,360° Rotatable	●	-	-
	Compensation Free Binocular,30°inclined,55mm-75mm	-	●	-
	Compensation FreeTrinocular,30°inclined,55mm-75mm	-	-	●
Eyeiece	WF10X/18mm	●	●	●
Objective	Achromatic 4X	●	●	●
	Achromatic 10X	●	●	●
	Achromatic 40X, Spring	●	●	●
	Achromatic 100X, Spring	●	●	●
Nosepiece	Quadruple Revolving Nosepiece	●	●	●
Stage	Double Layer Mechanical Stage,Size110mmx126mm,Moving Range:70mmx30mm	●	●	●
Condenser	ABBE N.A.1.25Condenser with Iris Diaphragm & Filter	●	●	●
Focusing	Coaxial Coarse & Fine Focusing Knobs,0.002mm,Focusing Range :30mm	●	●	●
Light Source	Halogen Bulb 6V/20W,,AC 220V/110V	●	●	●
Optional Accessories				
Head	Compensation Free Binocular,30°inclined,48mm-75mm			
	Compensation FreeTrinocular,30°inclined,48mm-75mm			
	Demonstration Head			
Eyeiece	WF16X			
	WF20X			
	P16X			
Objective	Achromatic 20X			
	Achromatic 60X(S)			
Light Source	LED			
	Reflective Mirror			
CCD Adapter	CCD Adapter			





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China