



1000X Laboratory Biological Microscope Dark Field Microscopes A11.1117

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1117
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month




Product Specification

- Applications: Laboratory
- Certification: CE
- Drawtube: Binocular, Trinocular
- Max Magnification: 1000X
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **high powered microscope ,
binocular biological microscope**



More Images



Product Description

Specifications:


A11.1117 Biological Microscope		-A	-B
Head	Compensation Free Binocular,30°inclined,55mm-75mm	●	-
	Compensation Free Trinocular,30°inclined,55mm-75mm	-	●
Eyepiece	WF10X/18mm	●	●
	Achromatic 4X	●	●
	Achromatic 10X	●	●

Our Product

Objective	Achromatic 40X, Spring	•	•
	Achromatic 100X, Spring,Oil	•	•
Nosepiece	Quadruple Revolving Nosepiece	•	•
Stage	Double Layer Mechanical Stage,Size150mmx150mm,Moving	•	•
	Range:80mmx50mm	•	•
Focusing	Coaxial Coarse and Fine Focusing Scale Value 0.002mm, Focusing Range :30mm	•	•
Condenser	ABBE N.A.1.25Condenser with Iris Diaphragm & Filter	•	•
Illumination	Kohler Illumination	•	•
Light Source	Halogen Bulb 6V/20W,AC 220V/110V,Brightness Adjustable	•	•
Optional Accessories			
Eyepiece	WF16X		
	WF20X		
	P16X		
Objective	Achromatic 20X		
	Achromatic 60X, Spring		
	Semi-plan Objective 4X		
	Semi-plan Objective 10X		
	Semi-plan Objective 20X		
	Semi-plan Objective 40X(S)		
	Semi-plan Objective 60X(S)		
	Semi-plan Objective 100X(S)		
	Plan Objective 4X		
	Plan Objective 10X		
	Plan Objective 20X		
	Plan Objective 40X(S)		
Plan Objective 60X(S)			
Plan Objective 100X(S)			
Dark Field	Darkfield		
Photo Adapter	Photo Adapter		
CCD Adapter	0.5x		
	0.57x		
	0.75x		



www.cnoec.com







Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China