



1000x Fine / Coarse Adjustment Microscope For School Laboratory A11.1121

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A11.1121
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School, Laboratory
- Certification: CE
- Drawtube: Binocular, Trinocular
- Max Magnification: 1000X
- Product Name: Microscope
- Theory: Biological Microscope
- Highlight: **binocular biological microscope,
biological compound microscope**



Product Description

1000x Fine / Coarse Adjustment Microscope For School Laboratory A11.1121

Specifications:

A11.1121 Biological Microscope		-A	-B
Head	Sliding Binocular, 45°inclined, 52mm-74mm	●	-
	Sliding Trinocular, 45°inclined, 52mm-74mm	-	●
Eyepiece	WF10X	●	●
Objective	Achromatic 4X	●	●
	Achromatic 10X	●	●
	Achromatic 40X, Spring	●	●
	Achromatic 100X, Spring, Oil	●	●
Nosepiece	Quadruple Revolving Nosepiece	●	●
Stage	Double Layer Mechanical Stage, Size 140x140mm, Moving Range: 75mmx50mm	●	●
Focusing	Coaxial Coarse and Fine Focusing Scale Value 0.002mm	●	●
Condenser	ABBE N.A. 1.25 Condenser with Iris Diaphragm & Filter	●	●
Collector	High Brightness Light Collector	●	●
Light Source	Halogen Bulb 6V/30W, AC 85V-230V, Adjustable Brightness	●	●
Optional Accessories			
Eyepiece	WF16X		
	WF20X		
	P16X		
Objective	Semi-plan Objective 4X		
	Semi-plan Objective 10X		
	Semi-plan Objective 20X		
	Semi-plan Objective 40X(S)		
	Semi-plan Objective 60X(S)		
	Semi-plan Objective 100X(S)		





www.cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China