



Monocular / Binocular Biological Microscope With Low Position Coaxial Coarse A12.0902

Our Product Introduction

for more products please visit us on cnoec.com

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.0902
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory, Hobby
- Certification: CE, Rohs
- Drawtube: Binocular, Monocular, Trinocular
- Illumination: Halogen Lamp
- Product Name: Compound Optical Microscope
- Theory: Biological Microscope
- Highlight: **binocular compound microscope, wide field microscope**



More Images



Product Description

Monocular / Binocular Biological Microscope With Low Position Coaxial Coarse A12.0902

Description:

- 195mm Achromatic Objectives, Mechanical Tube Length 160mm
- WF10x/18mm Large Field, Wide Angle Eyepiece, All Can Be Mounted With Reticle or Pointer
- Monocular, Binocular, Trinocular Head Optional
- Low Position Coaxial Coarse & Fine Focus System
- Double Layer Mechinal Stage

Specifications:

A12.0902 Laboratory Biological Microscope		A	B	C	D
Magnification	40x~1600x				
Head	Monocular Head, 45° Inclined, 360° Rotatable	o			
	Sliding Binocular Head, 45° Inclined, 360° Rotatable	o			
	Compensation Free Binocular Head, 30° Inclined, 360° Rotatable	o			
	Compensation Free Trinocular Head, 30° Inclined, 360° Rotatable	o			
Eyepiece	Wide Field WF10x/18mm, High Eyepoint	o	o	o	
	Wide Field WF16x/14mm, High Eyepoint				o
Nosepiece	Quadruple Nosepiece				
Objective	Achromatic 4x/0.10, 10x/0.25, 40xs/0.65	o			
	Achromatic 4x/0.10, 10x/0.25, 40xs/0.65, 100xs/1.25(Oil)		o		
	Plan Achromatic 4x/0.10, 10x/0.25, 40xs/0.65, 100xs/1.25(Oil)			o	o
Stage	Single Layer Plain Stage, 136*120mm, With Slide Clip	o			
	Double Layer Mechanical Stage, 132*140mm, Moving Range 50*76mm, Precision 0.1mm		o	o	o
Focusing	Low Position Coaxial Coarse & Fine Focusing. Fine Focusing Scale 0.002mm, With Tension Adjustment, With Up Limit				
Condenser	5 Holes Disc Condenser, N.A.0.65	o			
	N.A.1.25 Abbe Condenser, Iris Diaphragm, With Filter Holder		o	o	o
Filter	Blue, Yellow, Green				
Light Source	12V15W Halogen Light, 110V/220V				
A12.0902 Laboratory Biological Microscope Optional Accessories					
Eyepiece	WF10x/18mm with Reticle	A51.0905-1018R			
	WF10x/18mm with Pointer	A51.0906-1018			
	WF15x/13mm, High Eyepoint	A51.0902-1513			
	WF20x/10mm, High Eyepoint	A51.0902-2010			
Objective	Achromatic 60xs/0.85	A52.0901-60			
	Plan Achromatic 25xs/0.40	A52.0903-25			
Adapter	1.0x C-Mount	A55.0921-10			
	0.5x C-Mount, For 1/2" CCD, Focus Adjustable	A55.0921-05			
	0.35x C-Mount, For 1/3" CCD, Focus Adjustable	A55.0921-35			



A12.0902-A Laboratory Biological Microscope



A12.0902-B Laboratory Biological Microscope



A12.0902-C Laboratory Biological Microscope



A12.0902-D Laboratory Biological Microscope



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China