



Polarizing Compound Optical Microscope Fluorescent LED Microscopes A12.0906

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.0906
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School, Hobby, Home
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: LED
- Product Name: Compound Optical Microscope
- Theory: Biological Microscope
- Highlight: **binocular compound microscope, optical compound microscope**



for more products please visit us on cnoec.com

Product Description

Description:

- Infinity Optical System Optional
- PL10x/18mm, Wide Field, High Eyepoint
- Achromatic, Plan Achromatic, And Infinity Plan Objectives Optional
- Critical And Kohler Illumination Optional
- 90~240V Wide Range Voltage
- Dark Field, Polarizing, Phase Contrast Attachment Optional
- Fluorescent Attachment Optional, Upgrade to A16.0906 Fluorescent Microscope
- Upgrade to A31.0906 Digital Biological Microscope

Specifications:

A12.0906 Biological Microscope		A	B
Head	Binocular, 30° Inclined, 360° Rotatable, Interpupillary Range 54~75mm	o	
	Trinocular, 30° Inclined, 360° Rotatable, Interpupillary Range 54~75mm	o	
Eyepiece	PL10x/18mm, Wide Field, High Eyepoint		
Objective	Achromatic 4x/0.10, 10x/0.25, 40xs/0.65, 100xs/1.25 Oil		
Nosepiece	Reversed Quadruple		
Focusing	Coaxial Coarse & Fine Focusing, Low Position Design, Coarse Focusing Range 28mm, Fine Focusing Precision 0.002mm, With Tension Adjustment and Up Limit		
Working Stage	Double Layer Mechanical Stage, 150*140mm, Moving Range 76*50mm, Precision 0.1mm, Double Slide Holer, Low Position Coaxial Regulating Handles		
Condenser	Abbe N.A.1.25 Condenser, With Iris Diaphragm, With Filter Holder		
Filter	Blue, Dia. 45mm		
Light Source	Critical Illumination, 6V20W Halogen Lamp, Center Pre-adjusted		
Power	90~240V Wide Range Voltage		
A12.0906 Biological Microscope Optional Accessories			
Eyepiece	PL10x/20mm, Wide Field, High Eyepoint	A51.0903-1020	
	PL10x/18mm, Wide Field, High Eyepoint, With Pointer	A51.0905-1018R	
	PL10x/18mm, Wide Field, High Eyepoint, With Reseau	A51.0905-1018G	
	WF15x/13mm	A51.0902-1513	
	WF20x/10mm	A51.0902-2010	
	Plan Achromatic 4x/0.10	A52.0903-4	
	Plan Achromatic 10x/0.25	A52.0903-10	
	Plan Achromatic 20x/0.40	A52.0903-20	
	Plan Achromatic 40xs/0.65	A52.0903-40	

Objective	Plan Achromatic 100xs/1.25	A52.0903-100
	Infinity Plan Achromatic 4x/0.10	A52.0906-4
	Infinity Plan Achromatic 10x/0.25	A52.0906-10
	Infinity Plan Achromatic 20x/0.50	A52.0906-20
	Infinity Plan Achromatic 40x/0.65	A52.0906-40
	Infinity Plan Achromatic 100xs/1.25	A52.0906-100
Nosepiece	Reversed Quintuple Nosepiece	A54.0930-BS5
Filter	Yellow, Dia. 45mm	A56.0935-XSY
	Green, Dia. 45mm	A56.0935-XSG
	Netural, Dia. 45mm	A56.0935-XSF
Adapter	C-Mount 1.0x, Focus Adjustable	A55.0925-10
	C-Mount 0.5x , Focus Adjustable, For 1/2" CCD	A55.0925-05
	C-Mount 0.35x, Focus Adjustable , For 1/3" CCD	A55.0925-35
	Photo Adapter, With PK or MD Mount, With 3.2x Photo Eyepiece	A55.0901-2
Light Source	Kohler Illumination, 6V20W Halogen Lamp, Center Pre-adjusted	A56.0904-LK
	Critical Illumination, 3W LED, Center Pre-adjusted	A56.0904-LED
	Kohler Illumination, 3W LED, Center Pre-adjusted	A56.0904-LEDK
Condenser	Abbe N.A.1.25 Condenser, With Iris Diaphragm, With Filter Holder, With Polarizing & Phase Contrast Socket, For Critical Illumination	A56.0933-B
Condenser	Abbe N.A.1.25 Condenser, With Iris Diaphragm, With Filter Holder, With Polarizing & Phase Contrast Socket, For Kohler Illumination	A56.0933-BC
Dark Field	Dark Field Slide, For 4x~10x Objective	A5D.0915-1
	Dark Field Slide, For 20x~40x Objective	A5D.0915-2
Polarizing	Polarizer	A5P.0912-P
	Analyzer	A5P.0912-A
Plan Phase Contrast	Centering Telescope Eyepiece, Dia. 23.2mm	A5C.0910-232
	Phase Contrast Slide, For 10x, 40x Objective	A5C.0915-P1
	Phase Contrast Slide, For 20x, 100x Phase Contrast Objective	A5C.0915-P2
	Plan Phase Contrast Objective 10x/0.25	A5C.0931-10
	Plan Phase Contrast Objective 20x/0.50	A5C.0931-20
	Plan Phase Contrast Objective 40x/0.65	A5C.0931-40
	Plan Phase Contrast Objective 100x/1.25	A5C.0931-100

Infinity PlanPhase Contrast	Centering Telescope Eyepiece, Dia. 23.2mm	A5C.0910-232
	Phase Contrast Slide, For 10x, 40x Infinity Objective	A5C.0915-IP1
	Phase Contrast Slide, For 20x, 100x Infinity Objective	A5C.0915-IP2
	Infinity Plan Phase Contrast Objective 10x/0.25	A5C.0935-10
	Infinity Plan Phase Contrast Objective 20x/0.50	A5C.0935-20
	Infinity Plan Phase Contrast Objective 40x/0.65	A5C.0935-40
	Infinity Plan Phase Contrast Objective 100x/1.25	A5C.0935-100





A12.0906 Biological Microscope, Side View, Front View, Back View



Binocular Head, Diopter Adjustable



High Precision Achromatic Objectives



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China