



40X 100X Compound Optical Microscope Laboratory Biological Microscope A12.1004

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1004
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Student, School, Laboratory
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: Halogen Lamp, LED
- Product Name: Compound Optical Microscope
- Theory: Biological Microscope
- Highlight: **optical compound microscope, wide field microscope**



for more products please visit us on cnoec.com

Product Description

Compound Optical Microscope Laboratory Biological Microscope A12.1004

Specifications:


A12.1004 Biological Microscope		
Head	Seidentopf binocular viewing head, inclined at 30°, 360° rotatable	
Eyepiece	Wide Field Eyepiece WF10x/ 18	
Objective	Achromatic objective 4x	
	Achromatic objective 10x	
	Achromatic objective 40x	
	Achromatic objective 100x	
Nosepiece	Quadruple nosepiece	
Stage	Double layers mechanical stage, cross travel 78×54mm, using low-positioned X/ Y coaxial control knob	
Focusing	Coaxial Coarse & Fine knob	
Condenser	Abbe condenser N.A.1.25 with aperture diaphragm	
Illumination	3W LED	
Optional Accessories		Item No.
Head	Seidentopf trinocular viewing head, inclined at 30°, 360° rotatable	A53.1004-2
Objective	WF15x	A51.1002-15
	Plan achromatic objective 4x	A52.1004-4
	Plan achromatic objective 10x	A52.1004-10
	Plan achromatic objective 40x	A52.1004-40
	Plan achromatic objective 100x	A52.1004-100
	Infinite plan achromatic objective 4x	A52.1006-4
	Infinite plan achromatic objective 10x	A52.1006-10
	Infinite plan achromatic objective 40x	A52.1006-40
Photo Attachment		A55.1001
		A55.1002
Video Attachment		A55.1002
Dark Field Condensor		A5D.1001
Illumination	12V/ 20W Halogen lamp, brightness adjustable	A56.1003
	Kohler illumination	



www.optoec.com

 0086 13911110627

 sale@optoedu.com

 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China