



## 1000X Trinocular / Binocular Compound Microscope A12.1006

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1006
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Laboratory, Hobby
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: Halogen Lamp
- Max Magnification: 1000X
- Product Name: Compound Optical Microscope
- Highlight: **binocular compound microscope ,  
optical compound microscope**



## Product Description

### 1000X Trinocular / Binocular Compound Microscope A12.1006

#### Specifications:

A12.1006 Biological Microscope		A	B
Head	Siedentopf binocular head, inclined at 30°	●	
	Siedentopf trinocular head, inclined at 30°		●
Eyepiece	Wide field eyepiece WF10x/ 18	●	●
Objective	Semi-plan Achromatic Objective 4x	●	●
	Semi-plan Achromatic Objective 10x	●	●
	Semi-plan Achromatic Objective 40x	●	●
	Semi-plan Achromatic Objective 100x	●	●
Nosepiece	Backward quadruple nosepiece	●	●
Stage	Double Layers Mechanical Stage 140x140mm/ 75x50mm	●	●
Condenser	Abbe condenser NA.1.25 with Iris diaphragm & filter	●	●
Filter	Blue Filter	●	●
	Green Filter	●	●
Illumination	Halogen lamp 6V/ 20W, brightness adjustable	●	●
Optional Accessories		Item No.	
Eyepiece	Wide Field Eyepiece WF16x/ 11	A51.1002-16	
	Plan achromatic objective 4x	A52.1004-4	
Objective	Plan achromatic objective 10x	A52.1004-10	
	Plan achromatic objective 20x	A52.1004-20	
	Plan achromatic objective 40x	A52.1004-40	
	Plan achromatic objective 60x	A52.1004-60	
	Plan achromatic objective 100x	A52.1004-100	
Nosepiece	Quintuple nosepiece	A54.1001-5M	
Illumination	Kohler illumination	A56.1005-3	
Epi-fluorescent Attachment		A5F.1001	
Phase Contrast Attachment	Simple phase contrast set	A5C.1001	
	Sliding phase contrast set	A5C.1002-C1	
	Turret phase contrast set	A5C.1003-B1	
Polarization Set	Simple polarization set	A5P.1001	
Dark Field	Dark field condenser (Dry) N.A. 0.9	A5D.1001-BD	
Condenser	Dark field condenser (Oil) N.A. 1.25-1.36	A5D.1001-BO	






**Opto-Edu (Beijing) Co., Ltd.**

---

 0086 13911110627

 [sale@optoedu.com](mailto:sale@optoedu.com)

 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China