



Dark Field Compound Optical Microscope Trinocular With Blue / Green Filter CE A12.1008

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1008
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: School, Laboratory
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: Halogen Lamp
- Product Name: Compound Optical Microscope
- Theory: Biological Microscope
- Highlight: **binocular compound microscope , optical compound microscope**



for more products please visit us on cnoec.com

Product Description

Dark Field Compound Optical Microscope Trinocular With Blue / Green Filter CE A12.1008

Specifications:

A12.1008 Biological Microscope		A	B
Head	Siedentopf binocular head, inclined at 30°, 360° Rotatable	●	
	Siedentopf trinocular head, inclined at 30°, 360° Rotatable		●
Eyepiece	Wide field eyepiece WF10x/ 20		
Objective	Semi-plan Achromatic Objective 4x		
	Semi-plan Achromatic Objective 10x		
	Semi-plan Achromatic Objective 40x		
	Semi-plan Achromatic Objective 100x		
Nosepiece	Backward quadruple nosepiece		
Stage	Double Layers Mechanical Stage 140×140mm/ 75×50mm		
Condenser	Abbe condenser N.A.1.25 with Iris diaphragm & filter		
Filter	Blue filter		
	Green filter		
Illumination	Halogen lamp 6V/ 20W, brightness adjustable		
Optional Accessories		Item No.	
Eyepiece	Wide Field Eyepiece WF16x/ 11	A51.1002-16	
Objective	Plan Achromatic Objective 4x	A52.1004-4	
	Plan Achromatic Objective 10x	A52.1004-10	
	Plan Achromatic Objective 20x	A52.1004-20	
	Plan Achromatic Objective 40x	A52.1004-40	
	Plan Achromatic Objective 60x	A52.1004-60	
	Plan Achromatic Objective 100x	A52.1004-100	
Nosepiece	Quintuple nosepiece	A54.1001-5M	
Illumination	Kohler illumination	A56.1005-3	
Photo Attachment		A55.1001	
Video Attachment	C Mount 1x	A55.1001-10	
Epi-fluorescent Attachment	C Mount 0.5x	A55.1001-05	
Phase Contrast Attachment		A5F.1001	
Polarization Set	Simple phase contrast set	A5C.1001	
	Sliding phase contrast set	A5C.1002-C3	
	Turret phase contrast set	A5C.1003-B1	
Dark Field Condenser	Simple polarization set	A5P.1001	
Dark Field Condenser	Dark field condenser (Dry) N.A. 0.9	A5D.1001-BD	
	Dark field condenser (Immersion) N.A. 1.25-1.36	A5D.1001-BO	



www.cnoec.com



Easy Lamp Replacement Device



Holding Position In Back For Convenient Moving



Turret Phase Contrast Kit



Dry, Immersion Dark Field Condenser



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China