



Hobby Achromatic LED Wide Field Microscope Phase Contrast Light Microscopes A12.1010

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A12.1010
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory, Hobby
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: Halogen Lamp, LED
- Product Name: Compound Optical Microscope
- Theory: Biological Microscope
- Highlight: **optical compound microscope, wide field microscope**



for more products please visit us on cnoec.com

Product Description

Hobby Achromatic LED Wide Field Microscope Phase Contrast Light Microscopes A12.1010

Description:

- Universal Infinity System
- 3W LED High Light, or 6V20W Halogen

Specifications:

A12.1010 Biological Microscope		A	B
Head	Siedentopf binocular viewing head, inclined at 30°,360°rotatable	●	
	Siedentopf trinocular viewing head, inclined at 30°,360°rotatable		●
Optical System	Infinite optical system		
Eyepiece	WF10x/20		
Objective	Infinite plan achromatic objective 4x		
	Infinite plan achromatic objective 10x		
	Infinite plan achromatic objective 40x		
	Infinite plan achromatic objective 100x		
Nosepiece	Backward quadruple nosepiece		
Focusing	Coaxial coarse & fine adjustment, range 25mm, fine division 0.002mm		
Stage	Rectangular 216×150mm, cross travel 78×54mm, using low-positioned X/ Y coaxial control knob		
Condenser	Abbe condenser NA.1.25 with aperture diaphragm		
Illumination	12V/20W Halogen lamp, brightness adjustable		
Optional Accessories		Item No.	
Eyepiece	Wide Field Eyepiece WF15x/16	A51.1002-15B	
Nosepiece	Quintuple nosepiece	A54.1001-320M	
Illumination	3W LED	A56.1004-3	
Video	C-Mount adapter 1x	A55.1002-10	
Attachment	C-Mount adapter 0.5x	A55.1002-05	
Phase Contrast Attachment	Infinity plan phase contrast objective 10x/40x/100x+Condenser+Centering Telescope	A5C.1004	
Polarization Set		A5P.1004	





Seidentopf Trinocular Head, Inclined 30°, Rotatable 360°



Infinity Plan Achromatic Objectives 4x, 10x, 40x, 100x



Newly Designed X Direction Travel Guide To Avoid Scratch



Light Source Switch, Brightness Adjusting Knob



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China