



Halogen Lamp Inverted Optical Microscope Wide Field Microscopes CE A14.1101

Our Product Introduction


Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A14.1101
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Research
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: Halogen Lamp
- Product Name: Inverted Optical Microscope
- Theory: Inverted Microscope
- Highlight: **inverted fluorescence microscope, inverted light microscope**



More Images




for more products please visit us on cnoec.com

Product Description

Halogen Lamp Inverted Optical Microscope Wide Field Microscopes CE A14.1101

Specifications:

A14.1101 Inverted Metallurgical Microscope Specification		-A	-B	-C	-D
Head	Compensation Free Binocular,45°inclined,50mm-75mm	●	-	●	-
	Compensation Free Trinocular,45°inclined,50mm-75mm	-	●	-	●
Eyepiece	Wide field WF10X/20mm	●	●	-	-
	Wide field WF10X/23mm	-	-	●	●
Long Working Distance Plan-Objective	LWDPL 4X/0.1 W.D. 25mm	●	●	-	-
	LWDPL 10X/0.25 W.D. 11mm	●	●	-	-
	LWDPL 20X/0.4 W.D. 8mm	●	●	-	-
	LWDPL 40X/0.6 W.D. 3.8mm	●	●	-	-
Infinity Long Working Distance Objective	LWDPL 4X/0.1 W.D. 25mm	-	-	●	●
	LWDPL 10X/0.25 W.D. 12mm	-	-	●	●
	LWDPL 20X/0.4 W.D. 10mm	-	-	●	●
	LWDPL 40X/0.6 W.D. 7mm	-	-	●	●
Long Working Distance Phase Contrast Objective	LWDPL 20X	●	●	●	●
Focus System	Coaxial Coarse and Fine Focusing with Rack and Pinion Mechanism; Fine Focusing Scale Value 0.002mm.	●	●	●	●
Nosepiece	Quadruple	-	-	●	●
	Quintuple	●	●	-	-
Condenser	N.A.0.3 ABBE Condenser,W.D.75MM	●	●	●	●
Phase Contrast Sliding Plate	Center of Phase Contrast Ring Plate Adjustabel	●	●	●	●
Centering Telescope	Centering Telescope	●	●	●	●
Filter	Blue Filter	●	●	●	●
	Green Filter	●	●	●	●
	Yellow Filter	●	●	●	●
	Frosted Glass	●	●	●	●
Stage	Double Layer Mechanical, Central Stage: Φ110mm (Size242mmX172mm,Moving Range:75mmX50mm)	●	●	●	●
Light Source	Halogen Bulb,12V/30W,Adjustble Brightness	●	●	●	●
Optional Accessories					
Eyepiece	Wide field WF10X/18mm with Crosshair Graduation of 0.1mm Scale				
Long Working Distance Phase Contrast Objective	LWDPL 10X				
	LWDPL 40X				
Fluorescence Attachment	B Wave Band				
	G Wave Band				
	220V/110V/100W Mercury Lamp				
Photograph Attachment	Photography Connecting Tube,Adapter with MD or PK Mount,4X				
	Photography Eyepiece				





Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China