



Infinity Plan Fluorescence Microscope Halogen Lamp A16.1029 40x - 1000x

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A16.1029
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Application: Research, Learning
- Certification: CE, Rohs
- Drawtube: Trinocular
- Illumination: Halogen Lamp, LED
- Magnification: 40x - 1000x
- Theory: Fluorescence Microscope
- Highlight: **epi-fluorescence microscope, inverted fluorescence microscope**



Product Description

Infinity Plan Fluorescence Microscope Halogen Lamp A16.1029

Specifications:

A16.1029 Fluorescent Microscope					
Head	Seidentopf Trinocular Head, 30° inclined, interpupillary 48-75mm				
Eyepiece	EWF10X/22mm				
Nosepiece	Backward Quintuple Nosepiece				
Infinity Objective	Plan Achromatic 4X				
	Plan Achromatic 10X				
	Plan Achromatic 40X				
Stage	Plan Achromatic 100X				
	Double Layers Mechanical Stage 185x142mm, Moving Range 75x55mm				
Condenser	Swing Condenser N.A 0.9/0.25				
Focusing	Coaxial Coarse & Fine Focus Adjustment, Fine Division 0.001mm				
Illumination	External Kohler Illumination, Aspherical Collector, Halogen Lamp 6V30W				
Optional Accessories				Item No.	
Infinity Objective	Plan Achromatic 20X				
Attachement	Photo Attachement				
Video Adapter	Video Adapter				
Epi-fluorescent Attachment Specification					
Reflected Illumination	B(Blue)excitation filter system	Excitation	Dichroic Mirro	Barrier filter	
		BP460-490	DM505	BA515	
Lamp	3W LED(465-475nm)				
Immersion Oil	Fluorescence Free Oil				
Epi-fluorescent Attachment Optional Accessories					
Reflected Illumination	G(Green)Excitation Filter System	Excitation	Dichroic Mirro	Barrier filter	Item No.
		BP510 550	DM570	BA590	
Fluorescent Objective	Infinity Plan Fluorescent objective 4x				
	Infinity Plan Fluorescent objective 10x				
	Infinity Plan Fluorescent objective 20x				
	Infinity Plan Fluorescent objective 40x				
Lamp	Infinity Plan Fluorescent objective 100x				
	3W LED(520-530nm)				



 0086 13911110627  sale@optoedu.com  cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China