



## 1000x LED Fluorescence Microscope With 30°Inclined Head A16.1030

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A16.1030
- Minimum Order Quantity: 1 pc
- Price: Negotiation
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Application: Learning, Laboratory
- Certification: CE, Rohs
- Drawtube: Binocular, Trinocular
- Illumination: LED
- Magnification: 40x - 1000x
- Theory: Fluorescence Microscope
- Highlight: **inverted fluorescence microscope, epifluorescence microscope**



## Product Description

### 1000x LED Fluorescence Microscope With 30°Inclined Head A16.1030

#### Specifications:

A16.1030 Fluorescent Microscope				
		B	T	
Head	Seidentopf Binocular Head, 30° inclined, interpupillary 48-75mm	●	-	
	Seidentopf Trinocular Head, 30° inclined, interpupillary 48-75mm	-	●	
Eyepiece	WF10X/18mm			
Nosepiece	Backward Quadruple Nosepiece			
Infinity Objective	Semi-plan Achromatic 4X			
	Semi-plan Achromatic 10X			
	Semi-plan Achromatic 40X			
	Semi-plan Achromatic 100X			
Stage	Double Layers Mechanical Stage 140x140mm, Moving Range 75x50mm			
Condenser	Condenser N.A 1.25			
Focusing	Coaxial Coarse & Fine Focus Adjustment, Fine Division 0.002mm			
Illumination	LED Illumination, Brightness Adjustable.			
Optional Accessories			Item No.	
Infinity Objective	Semi-plan Achromatic 20X			
Attachment	Photo Attachment			
Video Adapter	Video Adapter			
Epi-fluorescent Attachment Specification				
Reflected Illumination	Excitation	Dichroic Mirro	Barrier filter	
B(Blue)excitation filter system	BP460-490	DM505	BA515	
Lamp	3W LED(465-475nm)			
Immersion Oil	Fluorescence Free Oil			
Epi-fluorescent Attachment Optional Accessories				
Reflected Illumination	Excitation	Dichroic Mirro	Barrier filter	Item No.
G(Green)Excitation Filter System	BP510 550	DM570	BA590	A5F.1010-G
Fluorescent Objective	Infinity Plan Fluorescent objective 4x			A5F.1031-4
	Infinity Plan Fluorescent objective 10x			A5F.1031-10
	Infinity Plan Fluorescent objective 20x			A5F.1031-20
	Infinity Plan Fluorescent objective 40x			A5F.1031-40
	Infinity Plan Fluorescent objective 100x			A5F.1031-100
Lamp	3W LED(520-530nm)			A56.1004-3




www.cnoec.com



**Opto-Edu (Beijing) Co., Ltd.**

 0086 13911110627

 [sale@optoedu.com](mailto:sale@optoedu.com)

 [cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China