



7x - 45x Medical Stereo Optical Microscope With 360°Rotatable Head A23.0902-S2

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.0902-S2
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Magnification: 7x - 45x
- Applications: Gem|Jewelry|Medical
- Certification: CE|Rohs
- Drawtube: Binocular|Trinocular
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,
stereo binocular microscopes**



for more products please visit us on cnoec.com

Product Description

7x - 45x Medical Stereo Optical Microscope With 360°Rotatable Head A23.0902-S2

Description:

- 0.7~4.5x Zoom Lens, with WF10x/20mm Eyepiece,
- Total Optical Magnification 7~45x
- Double Arm Universal Pole Stand

Specifications:

A23.0902-S2 Zoom Stereo Microscope Specification		B1	T1
Head	Binocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm	o	
	Trinocular Head, 45° inclined, 360°Rotatable, ±5 Diopter Adjustable, Interpupillary Adjustment 54-75mm	o	
Eyepiece	WF10x/22mm, High Eyepoint		
Zoom Lens	0.67~4.5x		
Zoom Ratio	1:6.7		
Magnification	6.7x~45x, (2x~225x with Auxliary Objective)		
Working Distance	100mm, With Auxliary Objectives 26~287mm		
Light Source	No Light Source		
Stand	Double Arm Universal Pole Stand STL2, Height 380mm, Focusing Pole Height 140mm, Level Move Range 235mm, Base 210*260*20mm,		

A23.0902-S2 Zoom Stereo Microscope Optical Data Sheet								
Auxiliary Lens	0.3x	0.4x	0.5x	0.75x	1.0x	1.5x	2.0x	
Working Distance	287mm	217mm	177mm	117mm	100mm	47mm	26mm	
WF 10x	Magnification	2.0-13.5x 18.0x	2.7-18.0x	3.4-22.5x	5.0-33.8x	6.7-45.0x	10.1-67.5x	13.4-90.0x
	View Field(mm)	109.5-16.3	82.1-12.2	65.7-9.8	43.8-6.5	32.8-4.9	21.9-3.3	16.4-2.4
WF 15x	Magnification	3.0-20.3x	4.0-27.0x	5.0-33.8x	7.5-50.6x	10.1-67.5x	15.1-101.3x	20.1-135x
	View Field(mm)	79.6-11.9	59.7-8.9	47.8-7.1	31.8-4.7	23.9-3.6	15.9-2.4	11.9-1.8
WF 20x	Magnification	4.0-27.0x	5.4-36.0x	6.7-45.0x	10.1-67.5x	13.4-90.0x	20.1-135.0x	26.8-180x
	View Field(mm)	59.7-8.9	44.8-6.7	35.8-5.3	12.9-2.4	17.9-2.7	11.9-1.8	9-1.3
WF 25x	Magnification	5.0-33.8x	6.7-45.0x	8.4-56.3x	12.6-84.4x	16.8-112.5x	25.1-168.8x	33.5-225x
	View Field(mm)	44.8-6.7	33.6-5.0	26.9-4.0	17.9-2.7	13.4-2.0	9.0-1.3	6.7-1.0

A23.0902-S2 Zoom Stereo Microscope Optional Accessories		Cata. No.
Head	Binocular Head + WF10x/20mm	A23.0902-BH
	Trinocular Head + WF10x/20mm	A23.0902-TH
Eyepiece	WF10x/20mm, High-eyepoint, Outside Dia. 30mm	A51.0923-10
	WF15x/15mm, High-eyepoint, Outside Dia. 30mm	A51.0923-15
	WF20x/10mm, High-eyepoint, Outside Dia. 30mm	A51.0923-20
	WF10x/20mm,with Reticle	A51.0933-10
	WF15x/15mm,with Reticle	A51.0933-15
	WF20x/10mm,with Reticle	A51.0933-20
Auxiliary Lens	0.5X,W.D. 165mm	A52.0928-05
	0.75X,W.D. 105mm	A52.0928-75
	1.5X,W.D. 45mm	A52.0928-15
	2.0X,W.D. 30mm	A52.0928-20
CCD Adapter	0.50x, For 1/2" CCD	A55.0933-1/2
	0.35x, For 1/3" CCD	A55.0933-1/3



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China