



Parallel Gem Stereo Optical Microscope Transmitted Light Microscopes Rohs A23.1004

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.1004
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Gem, Jewelry, Medical
- Certification: CE, Rohs
- Drawtube: Binocular
- Illumination: LED
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes,
stereo binocular microscopes**



More Images



for more products please visit us on cnoec.com

Product Description

Parallel Gem Stereo Optical Microscope Transmitted Light Microscopes Rohs A23.1004

Specifications:

A23.1004 Fluorescent Stereo Microscope		A	B	C
Optical System	Parallel optical zoom system	●	●	●
Head	Binocular head, 20° iclination	●	●	●
Eyepiece	Diopter adjustment eyepiece EW10X/22	●	●	●
Objective	Plan achromatic objective 1x	●	●	●
Zoom Objective	0.8x-5x			●
	0.8x-6.4x		●	
	0.8x-8x	●		
Zoom Ratio	1: 6			●
	1: 8		●	
	1: 10	●		
Fluorescent Attachment	GFP-B(EX480-500,DM505,BA510-560)	●	●	●
	GFP-L(EX460-500,DM505,BA510)			
Working Distance	78mm	●	●	●
Focusing	Focusing range 105mm	●	●	●
	Coarse focusing stand		●	●
	Coaxial coarse & fine focusing unit	●		
Illumination	Transmitted & reflected LED illumination, brightness adjustable	●	●	●
	100W Ultra hi-voltage spherical mercury lamp, power supplier NFP-1, digital display	●	●	●
Optional Accessories		Item No.		
Head	Tilting binocular head, 0°-35° Inclination	A53.1011-4		
Eyepiece	WF15x/16	A51.1027-1516		
	WF20x/12	A51.1027-2012		
	WF30x/18	A51.1027-30		
Objective	Plan Apo-chromatic objective 0.5x	A52.1023-05A		
	Plan achromatic objective 1x	A52.1023-05		
	Plan achromatic objective 2x	A52.1023-20		
	Achromatic objective 0.3x	A52.1023-03		
Fluorescent Attachment	GFP-B(EX430-450,DM455,BA465-495)			
	YEP-B(EX490-510,DM525,BA530-560)			
	G(EX515-550,DM5750,BA590)			
	RFP-B(EX525-555,DM570,BA585)			
Iris Diaphragm		A56.1021		
Stage	Mechanical Stage			
	Temperature Controlled Stage	A54.1004		
Stand	Universal Table Stand			
	D Stand			
Photo Attachment		A55.1001		
Video Attachment	CCD amount 0.5x	A55.1002-05		
	CCD amount 1.0x	A55.1002-10		







Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China