



Trinocular Stereo Optical Microscope With Halogen Lamp A23.1101

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A23.1101
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Industry, Jewelry
- Certification: CE, Rohs
- Drawtube: Trinocular
- Illumination: Halogen Lamp
- Product Name: Stereo Zoom Microscope
- Theory: Stereo Microscope
- Highlight: **stereo zoom microscopes, stereo binocular microscopes**



More Images



for more products please visit us on cnoec.com

Product Description

Trinocular Stereo Optical Microscope With Halogen Lamp A23.1101

Specifications:

A23.1101 Stereo Microscope Specification		-T	-AT	-BT	-CT
Head	Binocula,55-75mm	●	-	-	-
	Trinocular,55-75mm	-	●	●	●
Eyeiece	SWF10X/23mm	●	●	●	●
Objective	Zoom Range 0.63-5X,Zoom Ratio 1:8	●	●	●	●
Working Distance	110mm	●	●	●	●
Illumination	Top:Incident Light Halogen Bulb12V/15W,AC 85-265V	●	●	●	-
	Bottom:Transmission Light Halogen Bulb12V/15W12V/20W,AC 85-265V	●	●	-	-
	Ring Fluorescent Lamp	-	-	-	●
Optional Accessories					
Eyeiece	WF10X/25				
	WF15X/17				
	WF20X/12.5				
Auxiliary Objectiv	0.5X				
	1.5X				
	2X				
Lightsource	Ring Fluorescent Lamp				
Photo Adapter	Photo Adapter				
CCD	CCD Camera, Colour1/3" High Resolution 520 TV Line				



www.cnoec.com



A23.1101

A23.1101-T



A23.1101-AT

A23.1101-BT



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China