



USB 2.0 Microscope Accessories A59.4202

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.4202
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Microscope
- Certification: CE, Rohs
- Output: USB 2.0
- Product Name: Microscope Accessories
- Resolution: 5.0M
- Senser: CMOS
- Highlight: **compound microscope parts,
parts of the microscope**



for more products please visit us on cnoec.com

Product Description

Our Product Introduction

Specifications:

Specification	A59.4202-300	A59.4202-500
Output	USB 2.0	
Pixels	3.0M	5.0M
Senser	1/2" CMOS	1/2.5" CMOS
Image Resolution	2048*1536	2592*1944
Pixel Size	3.2umx3.2um	2.2umx2.2um
Frame Rate	5FPS	5FPS
SN Ratio	43db	40.5db
Dynamic Range	61db	66.5db
Sensitivity	1.0V/lux-sec	0.53V/lux-sec
Power	USB power supply	
Function	Auto/Manual Black Balance	
Driver & Software	Support DirectShow , Twain and LAbview. Including OCX, SDK and Measuring Software	
System	Windows 2000, WindowsXP,Windows 7,Windows Vista	
Size	50*50*46mm	
PCB Board	39*39mm	



www.cnoec.com

for more products please visit us on cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 choec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China