



Achromatic Phase Contrast Microscope A19.2702

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A19.2702
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Applications: Laboratory|Research
- Certification: CE|Rohs
- Drawtube: Binocular
- Illumination: LED
- Max Magnification: 40X - 1000X
- Theory: Phase Contrast Microscope
- Highlight: **phase contrast microscopes,
inverted phase contrast microscope**

Product Description

Our Product Introduction

Specifications:

Specification	A19.2702 Phase Contrast Microscope	B	T
Optical System	UCIS Infinity independent achromatic optical system		
Head	Seidentopf binocular head 30° inclined, 360 ° rotatable interpupillary distance 52-75mm	●	
	Seidentopf trinocular head 30° inclined, 360 ° rotatable interpupillary distance 52-75mm		●
Eyepiece	WF10X/20mm, high eye-point Centering Telescope		
Objective	Infinity Achromatic objective 4x		
	Infinity Achromatic phase-contrast objective 10x		
	Infinity Achromatic phase-contrast objective 40x(S)		
	Infinity Achromatic phase-contrast objective 100(S/Oil)		
Focus System	Coaxial coarse and fine focus mechanism with markings. Fine focus sensitivity 0.001mm		
Nosepiece	Inward quadruple nosepiece		
Condenser	N.A. 1.25 Phase contrast condenser, center adjustable		
Stage	Mechanical stage with double layers, 156mmx138mm, movement range 76mmx52mm		
Light source	6V 20W halogen illumination, brightness adjustable		
	6V 30W halogen illumination, brightness adjustable		
	3W LED, brightness adjustable		



for more products please visit us on cnoec.com



Opto-Edu (Beijing) Co., Ltd.

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China