



Infinity Achromatic LED Phase Contrast Microscope 40X - 1000X A19.2701

Our Product Introduction

Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A19.2701
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



Product Specification

- Certification: CE
- Drawtube: Binocular
- Illumination: LED
- Theory: Phase Contrast Microscope
- Highlight: **phase contrast light microscope,
phase contrast microscopes**

for more products please visit us on cnoec.com

Product Description

Infinity Achromatic LED Phase Contrast Microscope 40X - 1000X A19.2701

Specification	A19.2701 Phase Contrast Microscope	B	T
Optical System	UCIS Infinity independent achromatic optical system		
Head	Seidentopf binocular head 30° inclined, 360 ° rotatable interpupillary distance 52-75mm	•	
	Seidentopf trinocular head 30° inclined, 360 ° rotatable interpupillary distance 52-75mm		•
Eyepiece	WF10X/18mm, high eye-point		
Objective	Infinity Achromatic objective 4x		
	Infinity Achromatic phase-contrast objective 10x		
	Infinity Achromatic phase-contrast objective 40x(S)		
	Infinity Achromatic phase-contrast objective 100(S/Oil)		
Focus System	Coaxial coarse and fine focus mechanism with markings. Fine focus sensitivity 0.001mm		
Nosepiece	Quadruple nosepiece inward facing with positive click stops		
Condenser	N.A. 1.25 Abbe condenser with iris diaphragm, with phase-contrast slide		
Stage	Mechanical stage with double layers, 142mmx135mm, movement range 76mmx 52mm		
Phase Condenser Attachment	PH slide 10x, 40x, 100x		
Light source	6V 20W halogen illumination, brightness adjustable		
	3W LED, brightness adjustable		

