



## Binocular Jewelry Microscope With Zoom Ratio 1:4 A24.0401 For Research

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE,
- Model Number: A24.0401
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Certification: CE
- Drawtube: Binocular
- Highlight: **laboratory binocular microscope, gemological microscope**

for more products please visit us on [cnoec.com](http://cnoec.com)

## Product Description

### Binocular Jewelry Microscope With Zoom Ratio 1:4 A24.0401 For Research

A 24.0401 Specification of Gem Microscope	
Head	Binocular head, interpupillary distance: 55-75mm, diopter adjustment: $\pm 5$ diopter
Body	Optic head and stand tilt to a more natural 40° viewing angle
Eyepiece	WF10x
Zoom Ratio	1:4
Objective Zoom Range	10x-80x
Focusing	Focus range: 50mm
Working Distance	95mm
Illumination	Voltage Input 90-230V
	Adjustable iris diaphragm to control width of 2-42mm light beam
	Built-in optical-fiber illumination system Tungsten-halogen darkfield/brightfield illumination



Our Product Introduction

for more products please visit us on [cnoec.com](http://cnoec.com)



Opto-Edu (Beijing) Co., Ltd.



0086 13911110627



[sale@optoedu.com](mailto:sale@optoedu.com)



[cnoec.com](http://cnoec.com)

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China