



## 3.6 Inch LCD Digital Eyepiece Camera 5.0M CMOS 8x Digital Zoom A59.2301

Our Product Introduction

### Basic Information

- Place of Origin: China
- Brand Name: CNOEC, OPTO-EDU
- Certification: CE, Rohs
- Model Number: A59.2301
- Minimum Order Quantity: 1 pc
- Price: FOB \$1~1000, Depend on Order Quantity
- Packaging Details: Carton Packing, For Export Transportation
- Delivery Time: 5~20 Days
- Payment Terms: T/T, West Union, Paypal
- Supply Ability: 5000 pcs/ Month



### Product Specification

- Applications: Microscope
- Certification: CE
- Product Name: Microscope Accessories
- Resolution: 5.0M
- Sensor: CMOS
- Highlight: microscope accessory,  
compound microscope parts



for more products please visit us on cnoec.com

### More Images



## Product Description

### 3.6 Inch LCD Digital Eyepiece Camera 5.0M CMOS 8x Digital Zoom A59.2301

**A59.2301 3.6" LCD Digital Eyepiece Camera**, Designed for all kinds microscope, telescope use, sharing micro image with your friends and colleagues easily.

- 3.6" LCD TFT Screen
- 5.0M CMOS Digital Camera
- Output Image To SD Card, Mini SD Card, USB 2.0, AV

A59.2301 3.6" LCD Digital Eyepiece Camera	
LCD	3.6" TFT LCD Screen
Senser	2.0M CMOS, Interpolation To 5.0M
Resolution	2560*1920, 5.0M 2048*1536, 3.0M 1600*1200, 2.0M
Image Format	Jpeg, Avi
Memeory Slot	SD, Mini SD
Digital Zoom	8x
Data Output	SD, Mini SD, USB2.0, AV
Self-Timer	Auto Off At 10/15/20 Seconds
Power	DC 5V Power Adapter
Adapt	Ocular Dia.23.2mm For Biological Microscope Eyepiece Tube
Accessories	--Software CD, With User Manual --Ring Adaptors Dia.30mm For Stereo Microscope Eyepiece Tube --USB Cable --AV Cable --Power Adapter
Packing	Packing Box 12*14*18cm, G.W. 1.0 Kgs 30pcs/Ctn, Carton 45*39*62cm, G.W.16.5 Kgs





Packing Box, 12\*14\*18cm, G.W. 1.5 Kgs



Used with Microscope Eyepiece Tube



**Opto-Edu (Beijing) Co., Ltd.**

☎ 0086 13911110627

✉ sale@optoedu.com

🌐 cnoec.com

F-1501 Wanda Plaza, No. 18 Shijingshan Road, Beijing 100043, China